



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 9/8/2020

OVENTEMP IN Celsius(°C): 107
Time IN: 18:09
In Date: 09/04/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:15
Out Date: 09/05/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB110845

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
L3856-01	MEW267	1	1.17	9.7	9.48	97.4	
L3856-02	MEW268	2	1.16	9.98	9.8	98.0	
L3856-03	MEW269	3	1.17	9.77	9.69	99.1	
L3856-04	MEW269D	4	1.17	9.77	9.69	99.1	
L3856-05	MEW269S	5	1.17	9.77	9.69	99.1	
L3856-06	MEW270	6	1.17	9.99	9.79	97.7	
L3856-07	MEW271	7	1.18	9.77	9.61	98.1	
L3856-08	MEW272	8	1.17	9.81	9.59	97.5	
L3856-09	MEW273	9	1.17	9.77	9.62	98.3	
L3856-10	MEW274	10	1.16	9.94	9.73	97.6	
L3856-11	MEW275	11	1.16	9.73	8.89	90.2	
L3856-12	MEW276	12	1.15	9.99	9.61	95.7	
L3856-13	MEW277	13	1.16	9.91	7.85	76.5	
L3856-14	MEW278	14	1.16	9.93	9.66	96.9	
L3856-15	MEW284	15	1.17	9.8	9.45	95.9	
L3856-16	MEW279	16	1.17	9.91	9.21	92.0	
L3856-17	MEW280	17	1.18	9.96	8.48	83.1	
L3856-18	MEW283	18	1.18	9.84	9.51	96.2	

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$