



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 9/10/2020

OVENTEMP IN Celsius(°C): 108
Time IN: 15:00
In Date: 09/09/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 105
Time OUT: 09:40
Out Date: 09/10/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB110879

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
L3945-02	BFWX6	1	1.14	9.76	8.45	84.8	
L3945-03	BFWX7	2	1.18	9.76	8.47	85.0	
L3945-04	BFWX8	3	1.15	9.96	8.92	88.2	
L3945-05	BFWX9	4	1.14	9.98	8.79	86.5	
L3945-06	BFWY0	5	1.18	9.57	8.23	84.0	
L3945-07	BFWY1	6	1.13	9.65	8.61	87.8	
L3945-08	BFWY2	7	1.1	9.71	8.97	91.4	
L3945-09	BFWY3	8	1.14	9.67	8.24	83.2	
L3945-10	BFWY4	9	1.13	9.61	8.67	88.9	
L3945-11	BFX00	10	1.18	9.94	8.71	86.0	
L3945-12	BFX04	11	1.12	9.92	9.67	97.2	
L3945-13	BFX05	12	1.16	9.69	8.83	89.9	
L3945-14	BFX06	13	1.1	9.65	9.24	95.2	
L3945-15	BFX07	14	1.17	9.62	8.3	84.4	
L3945-16	VHBLK01	15	1.00	2.00	2.00	100.0	VHBLK

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$