



PERCENT SOLID

Supervisor: apatel
Analyst: JIGNESH
Date: 9/16/2020

OVENTEMP IN Celsius(°C): 107
Time IN: 14:25
In Date: 09/15/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:05
Out Date: 09/16/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB110970

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
L3968-01	MYAQM6	1	1.15	9.9	8.1	79.4	
L3968-02	MYAQM5	2	1.18	9.61	7.4	73.8	
L3968-03	MYAQM4	3	1.15	9.95	7.88	76.5	
L3968-04	MYAQM3	4	1.13	9.96	8.54	83.9	
L3968-05	MYAQM2	5	1.14	9.93	7.87	76.6	
L3968-06	MYAQN5	6	1.18	9.8	8.16	81.0	
L3968-07	MYAQN4	7	1.17	9.96	6.8	64.1	
L3968-08	MYAQN3	8	1.15	9.98	8.17	79.5	
L3968-09	MYAQN2	9	1.12	9.59	7.69	77.6	
L3968-10	MYAQN1	10	1.18	9.89	7.85	76.6	
L3968-11	MYAQN0	11	1.13	9.7	8.55	86.6	
L3968-12	MYAQM9	12	1.15	9.96	6.92	65.5	
L3968-13	MYAQM8	13	1.1	9.91	7.97	78.0	
L3968-14	MYAQM7	14	1.14	9.97	8.12	79.0	
L3968-15	MYAQS0	15	1.14	9.96	4.12	33.8	
L3968-16	MYAQS1	16	1.16	9.95	7.08	67.3	
L3968-17	MYAQM1	17	1.13	9.64	7.85	79.0	
L3968-18	MYAQM1D	18	1.13	9.64	7.85	79.0	
L3968-19	MYAQM1S	19	1.13	9.64	7.85	79.0	

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$