



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 9/18/2020

OVENTEMP IN Celsius(°C): 107
Time IN: 17:25
In Date: 09/17/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:22
Out Date: 09/18/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB111029

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
L3870-06	SU-02-091720	1	1.17	9.67	9.22	94.7	
L3870-07	SU-02-091720	2	1.16	9.97	9.28	92.2	
L3870-08	SU-02-091720-1	3	1.1	9.51	9.00	93.9	
L3870-09	SU-02-091720-2	4	1.15	9.53	9.03	94.0	
L3871-13	TR-05-091720	5	1.13	9.92	8.85	87.8	
L3871-14	TR-05-091720	6	1.15	9.74	8.36	83.9	
L3871-15	TR-05-091720-1	7	1.16	9.97	8.97	88.6	
L3871-19	TR-08-091720	8	1.16	9.62	8.71	89.2	
L3871-20	TR-08-091720	9	1.13	9.83	8.8	88.2	
L3871-21	TR-08-091720-1	10	1.12	9.88	9.00	90.0	
L3930-06	VNJ-242	11	1.11	9.6	7.07	70.2	
L4017-01	SS-1	12	1.13	9.97	9.2	91.3	
L4017-02	SS-2	13	1.16	9.61	8.69	89.1	
L4017-03	SS-3	14	1.13	9.51	8.18	84.1	
L4017-04	SS-4	15	1.16	9.85	8.98	90.0	
L4017-05	SS-5	16	1.17	9.63	8.42	85.7	
L4017-06	SS-6	17	1.15	9.85	9.12	91.6	
L4017-07	SS-7	18	1.13	9.82	9.21	93.0	
L4017-08	SS-8	19	1.16	9.77	8.94	90.4	
L4017-09	SS-9	20	1.16	9.56	8.7	89.8	
L4017-10	SS-10	21	1.17	9.88	8.75	87.0	
L4017-11	SS-11	22	1.17	9.89	8.64	85.7	
L4034-06	SOIL-SP-MH-C01-MH-J09-02	23	1.13	9.62	8.83	90.7	
L4051-01	KTS-SB-02-1-2	24	1.13	9.85	8.71	86.9	
L4051-02	KTS-SB-01-1-2	25	1.16	9.87	8.1	79.7	
L4051-03	KTS-IB-01-1-2	26	1.11	9.97	8.9	87.9	
L4051-04	KTS-IB-02-1-2	27	1.16	9.85	8.92	89.3	
L4051-05	KTS-ST-02-1-2	28	1.13	9.88	9.36	94.1	



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 9/18/2020

OVENTEMP IN Celsius(°C): 107
Time IN: 17:25
In Date: 09/17/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:22
Out Date: 09/18/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB111029

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
L4051-06	KTS-ST-01-1-2	29	1.16	9.86	8.95	89.5	
L4052-01	HEF317G-TOP	30	1.00	2.00	2.00	100.0	CABLE SAMPLE
L4052-02	HED160B-TOP	31	1.00	2.00	2.00	100.0	CABLE SAMPLE
L4052-03	HEF317G-BOTTOM	32	1.00	2.00	2.00	100.0	CABLE SAMPLE
L4052-04	HED160B-BOTTOM	33	1.00	2.00	2.00	100.0	CABLE SAMPLE
L4052-06	CO-76728	34	1.00	2.00	2.00	100.0	OIL SAMPLE
L4053-01	1813	35	1.00	2.00	2.00	100.0	OIL SAMPLE
L4053-03	VNJ-231	36	1.00	2.00	2.00	100.0	OILY DEBRIS
L4053-04	VNJ-231	37	1.00	2.00	2.00	100.0	OILY DEBRIS
L4054-01	KZZ075N-TOP	38	1.00	2.00	2.00	100.0	CABLE SAMPLE
L4054-02	KZZ075N-BOTTOM	39	1.00	2.00	2.00	100.0	CABLE SAMPLE
L4057-01	JA-DRUM-COMP	40	1.00	2.00	2.00	100.0	OILY DEBRIS
L4057-02	JA-DRUM-COMP	41	1.00	2.00	2.00	100.0	OILY DEBRIS
L4059-01	PAV-B33-091720-0-1	42	1.16	9.6	9.06	93.6	
L4059-03	PAV-B33-091720-0-10	43	1.15	9.97	9.3	92.4	
L4059-05	PAV-B33-091720-10-19	44	1.17	9.93	9.19	91.6	
L4059-07	PAV-B33-091720-18-19	45	1.16	9.94	9.16	91.1	
L4062-01	EROC-CF-1	46	1.17	9.5	8.85	92.2	
L4062-02	EROC-CF-2	47	1.16	9.73	9.24	94.3	
L4062-03	EROC-NI4-1	48	1.17	9.96	9.83	98.5	
L4062-04	EROC-NI4-2	49	1.16	9.64	9.42	97.4	
L4067-01	PCB-091620-SR-01	50	1.00	2.00	2.00	100.0	CAULKING SAMPLE
L4067-02	PCB-091620-SR-02	51	1.00	2.00	2.00	100.0	CAULKING SAMPLE
L4068-01	PE-1A- (1-1.5)	52	1.18	9.67	8.68	88.3	
L4068-02	PE-2A- (1-1.5)	53	1.14	9.63	8.69	88.9	
L4068-03	PE-3A- (1-1.5)	54	1.16	9.82	9.04	91.0	
L4068-04	PE-4A- (1-1.5)	55	1.16	9.89	8.89	88.5	
L4068-05	PE-5- (5-5.5)	56	1.17	9.81	8.17	81.0	



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 9/18/2020

OVENTEMP IN Celsius(°C): 107
Time IN: 17:25
In Date: 09/17/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:22
Out Date: 09/18/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB111029

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
L4068-06	PE-1B-(4.5-5)	57	1.17	9.97	8.73	85.9	
L4068-07	PE-2B-(4.5-5)	58	1.16	9.98	8.89	87.6	
L4068-08	PE-3B-(4.5-5)	59	1.18	9.55	8.92	92.5	
L4068-09	PE-4B-(4.5-5)	60	1.16	9.96	8.71	85.8	
L4068-10	DUPLICATE	61	1.17	9.86	8.9	89.0	
L4073-01	H-1	62	1.18	9.8	9.26	93.7	
L4073-02	H-1A	63	1.17	9.55	9.00	93.4	
L4074-01	KTS-IB-0102-0-6	64	1.18	9.98	7.83	75.6	
L4074-02	KTS-ST-0102-0-6	65	1.17	9.85	8.35	82.7	
L4075-01	KTS-WC-S-IB-ST	66	1.18	9.83	7.93	78.0	

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$