



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 9/21/2020

OVENTEMP IN Celsius(°C): 107
Time IN: 14:25
In Date: 09/18/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:00
Out Date: 09/19/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB111040

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
L4002-01	MJLTQ0	1	1.13	9.66	8.02	80.8	
L4002-02	MJLTQ2	2	1.16	9.9	8.45	83.4	
L4002-03	MJLTQ4	3	1.16	9.53	8.69	90.0	
L4002-04	MJLTQ5	4	1.12	9.88	7.97	78.2	
L4002-05	MJLTQ8	5	1.18	9.79	7.3	71.1	
L4002-06	MJLTQ9	6	1.18	9.86	7.97	78.2	
L4002-07	MJLTR0	7	1.11	9.66	8.2	82.9	
L4002-08	MJLTR1	8	1.15	9.71	8.00	80.0	
L4002-09	MJLTR2	9	1.18	9.74	8.6	86.7	
L4002-10	MJLTR3	10	1.1	9.84	8.58	85.6	
L4002-11	MJLTR4	11	1.13	9.63	7.86	79.2	
L4002-12	MJLTR5	12	1.18	9.69	8.62	87.4	
L4002-13	MJLTT0	13	1.18	9.56	7.58	76.4	
L4002-14	MJLTT1	14	1.18	9.68	7.73	77.1	
L4002-15	MJLTT6	15	1.1	9.64	8.43	85.8	
L4002-16	MJLTT7	16	1.18	9.86	8.2	80.9	
L4002-17	MJLTY3	17	1.13	9.77	7.38	72.3	
L4002-18	MJLTY4	18	1.13	9.78	8.06	80.1	
L4002-19	MJLTX8	19	1.16	9.56	6.06	58.3	
L4002-20	MJLTX8D	20	1.16	9.56	6.06	58.3	
L4002-21	MJLTX8S	21	1.16	9.56	6.06	58.3	

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$