



PERCENT SOLID

Supervisor: apatel
Analyst: JIGNESH
Date: 9/22/2020

OVENTEMP IN Celsius(°C): 107
Time IN: 17:00
In Date: 09/21/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:22
Out Date: 09/22/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB111071

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
L3864-08	CL-02-092120	43	1.15	9.95	9.28	92.4	
L3864-09	CL-02-092120	44	1.16	9.77	9.22	93.6	
L3864-10	CL-02-092120-1	45	1.18	9.67	9.25	95.1	
L3864-11	CL-02-092120-2	46	1.13	9.62	9.1	93.9	
L3864-12	CL-02-092120-3	47	1.16	9.62	9.09	93.7	
L3864-13	CL-02-092120-4	48	1.13	9.59	8.61	88.4	
L3864-14	CL-02-092120-5	49	1.13	9.61	8.85	91.0	
L3866-06	OK-01-092120	59	1.13	9.73	9.46	96.9	
L3866-07	OK-01-092120	50	1.13	9.59	9.35	97.2	
L3866-08	OK-01-092120-1	51	1.13	9.83	9.28	93.7	
L3866-09	OK-01-092120-2	52	1.13	9.91	9.39	94.1	
L3872-07	BU-02-092120	53	1.14	9.98	9.33	92.6	
L3872-08	BU-02-092120	54	1.13	9.86	9.28	93.4	
L3872-09	BU-02-092120-1	55	1.13	9.68	9.00	92.0	
L4098-01	31232	56	1.1	9.93	9.82	98.8	
L4099-01	029	29	1.16	9.98	3.00	20.9	
L4099-02	029	30	1.18	9.52	3.35	26.0	
L4099-03	036	31	1.15	9.94	4.88	42.4	
L4099-04	036	32	1.1	9.73	3.79	31.2	
L4099-05	031	33	1.13	9.75	6.47	61.9	
L4099-06	031	34	1.13	9.54	5.06	46.7	
L4099-07	105	35	1.13	9.66	6.72	65.5	
L4099-08	105	36	1.15	9.96	5.61	50.6	
L4100-01	024	37	1.16	9.8	8.61	86.2	
L4100-02	024	38	1.13	9.68	8.47	85.8	
L4100-03	RBR-59	39	1.11	9.67	5.35	49.5	
L4100-04	RBR-59	40	1.16	9.97	5.6	50.4	
L4100-05	021	41	1.16	9.98	9.1	90.0	



PERCENT SOLID

Supervisor: apatel
Analyst: JIGNESH
Date: 9/22/2020

OVENTEMP IN Celsius(°C): 107
Time IN: 17:00
In Date: 09/21/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:22
Out Date: 09/22/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB111071

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
L4100-06	021	42	1.13	9.88	6.39	60.1	
L4101-01	0-6FT	57	1.13	9.86	9.33	93.9	
L4101-04	2-4FT	58	1.19	9.51	9.2	96.3	
L4102-04	T-1	1	1.00	2.00	2.00	100.0	OIL SAMPLE
L4102-05	PE-1A- (0.0-0.5)	2	1.15	9.77	7.45	73.1	
L4102-06	PE-1B- (3.0-3.5)	3	1.17	9.7	8.98	91.6	
L4102-07	PE-2A- (0.5-1.0)	4	1.13	9.55	8.26	84.7	
L4102-08	PE-2B- (3.0-3.5)	5	1.14	9.57	9.03	93.6	
L4102-09	PE-3A- (0.5-1.0)	6	1.12	9.72	8.93	90.8	
L4102-10	PE-3B- (3.0-3.5)	7	1.19	9.81	9.58	97.3	
L4102-11	PE-4A- (0.5-1.0)	8	1.14	9.81	9.06	91.3	
L4102-12	PE-4B- (3.0-3.5)	9	1.14	9.98	9.7	96.8	
L4102-13	PE-5- (3.5-4.0)	10	1.17	9.65	9.39	96.9	
L4102-14	DUP-092020	11	1.15	9.94	9.36	93.4	
L4103-01	PE-1- (1.5-2)	12	1.13	9.72	8.63	87.3	
L4103-02	PE-2- (2-2.5)	13	1.13	9.69	9.02	92.2	
L4103-03	PE-3- (1.0-1.5)	14	1.17	9.67	8.6	87.4	
L4103-04	PE-4- (0.5-1)	15	1.13	9.82	8.63	86.3	
L4103-05	PE-5- (2.5-3)	16	1.1	9.72	8.59	86.9	
L4103-06	PE-6- (2.5-3)	17	1.13	9.71	8.44	85.2	
L4103-07	PE-7- (2.5-3)	18	1.16	9.97	8.43	82.5	
L4103-08	PE-8- (2.5-3)	19	1.12	9.62	8.93	91.9	
L4103-09	PE-9- (2.5-3)	20	1.17	9.75	8.67	87.4	
L4103-10	DUPLICATE	21	1.19	9.85	8.36	82.8	
L4103-11	T-1	22	1.00	2.00	2.00	100.0	OIL SAMPLE
L4104-01	PAV-B29-092110-0-1	23	1.15	9.5	9.2	96.4	
L4104-03	PAV-B29-092110-0-10	24	1.14	9.54	8.72	90.2	
L4104-05	PAV-092110-DP	25	1.16	9.91	9.51	95.4	



PERCENT SOLID

Supervisor: apatel
Analyst: JIGNESH
Date: 9/22/2020

OVENTEMP IN Celsius(°C): 107
Time IN: 17:00
In Date: 09/21/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:22
Out Date: 09/22/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB111071

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
L4104-07	PAV-092110-DP	26	1.16	9.76	8.9	90.0	
L4104-09	PAV-B29-092110-10-17	27	1.19	9.86	8.89	88.8	
L4104-11	PAV-B29-092110-16-17	28	1.13	9.91	9.02	89.9	

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$