



PERCENT SOLID

Supervisor: apatel
Analyst: JIGNESH
Date: 9/28/2020

OVENTEMP IN Celsius(°C): 107
Time IN: 15:03
In Date: 09/25/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:03
Out Date: 09/26/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB111162

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
L4156-13	X2P17	1	1.00	2.00	2.00	100.0	P.E. SAMPLE
L4156-14	X2P17D	2	1.00	2.00	2.00	100.0	P.E. SAMPLE
L4156-15	X2P17S	3	1.00	2.00	2.00	100.0	P.E. SAMPLE
L4156-16	X2P18	4	1.00	2.00	2.00	100.0	P.E. SAMPLE
L4156-17	X2P18D	5	1.00	2.00	2.00	100.0	P.E. SAMPLE
L4156-18	X2P18S	6	1.00	2.00	2.00	100.0	P.E. SAMPLE
L4156-19	X2P19	7	1.00	2.00	2.00	100.0	P.E. SAMPLE
L4156-20	X2P19D	8	1.00	2.00	2.00	100.0	P.E. SAMPLE
L4156-21	X2P19S	9	1.00	2.00	2.00	100.0	P.E. SAMPLE
L4156-22	X2P20	10	1.00	2.00	2.00	100.0	P.E. SAMPLE
L4156-23	X2P20D	11	1.00	2.00	2.00	100.0	P.E. SAMPLE
L4156-24	X2P20S	12	1.00	2.00	2.00	100.0	P.E. SAMPLE

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$