



PERCENT SOLID

Supervisor: apatel
Analyst: JIGNESH
Date: 9/29/2020

OVENTEMP IN Celsius(°C): 107
Time IN: 17:00
In Date: 09/28/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:00
Out Date: 09/29/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB111191

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
L3867-07	EO-01-092820	34	1.13	9.76	9.08	92.1	
L3867-08	EO-01-092820	35	1.12	9.64	9.33	96.4	
L3867-09	EO-01-092820-1	36	1.13	9.56	9.00	93.4	
L3867-10	EO-01-092820-2	37	1.14	9.64	9.19	94.7	
L3867-11	EO-01-092820-3	38	1.13	9.91	9.08	90.5	
L3869-11	OR-03-092820	39	1.11	9.83	9.28	93.7	
L3869-12	OR-03-092820	40	1.14	9.77	8.83	89.1	
L3869-13	OR-03-092820-1	41	1.16	9.75	9.09	92.3	
L3869-14	OR-03-092820-2	42	1.13	9.59	9.12	94.4	
L4157-01	KMA9581-END-1	1	1.00	2.00	2.00	100.0	CABLE SAMPLE
L4157-02	KMA9581-END-2	2	1.00	2.00	2.00	100.0	CABLE SAMPLE
L4161-01	KMH490H-END-1	3	1.00	2.00	2.00	100.0	CABLE SAMPLE
L4161-02	KMH490H-END-2	4	1.00	2.00	2.00	100.0	CABLE SAMPLE
L4161-03	JEH607V-END-1	5	1.00	2.00	2.00	100.0	CABLE SAMPLE
L4161-04	JEH607V-END-2	6	1.00	2.00	2.00	100.0	CABLE SAMPLE
L4161-05	KMA696E-END-1-1-1	7	1.00	2.00	2.00	100.0	CABLE SAMPLE
L4161-06	KMA696E-END-2-1-2	8	1.00	2.00	2.00	100.0	CABLE SAMPLE
L4161-07	KMA696E-END-1-2-1	9	1.00	2.00	2.00	100.0	CABLE SAMPLE
L4161-08	KMA696E-END-2-2-2	10	1.00	2.00	2.00	100.0	CABLE SAMPLE
L4161-09	JEH735C-END-1	11	1.00	2.00	2.00	100.0	CABLE SAMPLE
L4161-10	JEH735C-END-2	12	1.00	2.00	2.00	100.0	CABLE SAMPLE
L4170-01	RBR251606	13	1.13	9.62	6.54	63.7	
L4170-02	RBR251606	14	1.15	9.65	8.48	86.2	
L4172-01	B-6(3-4)	15	1.12	9.83	8.73	87.4	
L4172-02	B-6(7-9)	16	1.16	9.97	8.91	88.0	
L4172-04	B-6(11-13)	17	1.16	9.6	8.37	85.4	
L4172-05	DUP-2	18	1.14	9.94	8.91	88.3	
L4172-07	B-7(0-1)	19	1.15	9.91	9.52	95.5	



PERCENT SOLID

Supervisor: apatel
Analyst: JIGNESH
Date: 9/29/2020

OVENTEMP IN Celsius(°C): 107
Time IN: 17:00
In Date: 09/28/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:00
Out Date: 09/29/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB111191

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
L4172-09	B-7 (9-11)	20	1.12	9.79	8.13	80.9	
L4172-10	B-7 (12-14)	21	1.15	9.96	8.57	84.2	
L4172-12	B-8 (9-11)	22	1.11	9.54	8.61	89.0	
L4172-13	B-8 (11-13)	23	1.16	9.7	8.69	88.2	
L4172-14	B-8 (13-14)	24	1.13	9.82	8.69	87.0	
L4173-01	TP-1- (3FT)	25	1.12	7.81	6.73	83.9	
L4173-02	TP-1- (0-6FT)	26	1.13	9.75	8.39	84.2	
L4173-04	TP-1- (0-6FT) -EPH-2	27	1.16	9.94	8.64	85.2	
L4173-05	TP-2- (3FT)	28	1.17	9.95	9.02	89.4	
L4173-06	TP-2- (0-6FT)	29	1.18	9.77	8.96	90.6	
L4173-08	TP-2- (0-6FT) -EPH-2	30	1.13	9.82	9.09	91.6	
L4173-09	TP-3- (3FT)	31	1.16	9.69	8.96	91.4	
L4173-10	TP-3- (0-6FT)	32	1.13	9.7	8.82	89.7	
L4173-12	TP-3- (0-6FT) -EPH-2	33	1.13	9.53	8.66	89.6	

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$