



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 9/30/2020

OVENTEMP IN Celsius(°C): 107
Time IN: 17:00
In Date: 09/29/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:10
Out Date: 09/30/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB111210

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
L4176-01	AC-4	1	1.12	9.55	9.1	94.7	
L4179-01	B-9 (2-3)	2	1.17	9.73	9.1	92.6	
L4179-02	B-9 (7-9)	3	1.13	9.84	9.11	91.6	
L4179-03	B-9 (12-14)	4	1.11	9.55	8.64	89.2	
L4179-04	B-10 (1-2)	5	1.15	9.89	9.5	95.5	
L4179-05	B-10 (4-5)	6	1.18	9.88	9.13	91.4	
L4179-06	B-10 (12-14)	7	1.14	9.82	8.6	85.9	
L4179-07	B-11 (0-1)	8	1.18	9.81	9.5	96.4	
L4179-08	B-11 (1-2)	9	1.18	9.93	9.6	96.2	
L4179-09	L4179-08MS	10	1.18	9.93	9.6	96.2	
L4179-10	L4179-08MSD	11	1.18	9.93	9.6	96.2	
L4179-11	B-11 (12-14)	12	1.13	9.93	8.86	87.8	
L4179-12	B-11 (9-11)	13	1.11	9.89	8.83	87.9	
L4189-01	TP-C-2334-C	14	1.12	9.74	8.7	87.9	

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$