



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 10/1/2020

OVENTEMP IN Celsius(°C): 107
Time IN: 17:00
In Date: 09/30/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:15
Out Date: 10/01/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB111222

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
L3863-09	HD-02-093020	42	1.14	9.96	8.47	83.1	
L3863-10	HD-02-093020	43	1.18	9.96	9.05	89.6	
L3863-11	HD-02-093020-1	44	1.18	9.95	8.67	85.4	
L3870-16	SU-03-093020	45	1.18	9.97	8.9	87.8	
L3870-17	SU-03-093020	46	1.14	9.97	9.6	95.8	
L3870-18	SU-03-093020-1	47	1.14	9.76	8.02	79.8	
L3870-19	SU-03-093020-2	48	1.12	9.66	8.59	87.5	
L4190-01	PAV-B14-092920-0-1	1	1.15	9.96	9.66	96.6	
L4190-03	PAV-B14-092920-0-10	2	1.12	9.53	8.99	93.6	
L4190-05	PAV-B14-092920-10-16	3	1.12	9.95	9.1	90.4	
L4190-07	PAV-B14-092920-15-16	4	1.13	9.85	8.66	86.4	
L4190-09	PAV-B18-092920-0-1	5	1.13	9.72	8.96	91.2	
L4190-11	PAV-B18-092920-0-10	6	1.15	9.98	9.33	92.6	
L4190-13	PAV-B18-092920-10-20	7	1.17	9.98	9.24	91.6	
L4190-15	PAV-B18-092920-19-20	8	1.13	9.9	8.42	83.1	
L4190-17	PAV-B18-092920-DP	9	1.11	9.8	9.06	91.5	
L4190-19	PAV-B18-092920-DP	10	1.16	9.66	8.95	91.6	
L4191-01	SS-34	11	1.14	9.64	8.99	92.4	
L4191-02	SS-33	12	1.13	9.58	8.91	92.1	
L4191-03	SS-83	13	1.15	9.94	9.26	92.3	
L4191-04	SS-32	14	1.13	9.8	9.4	95.4	
L4191-05	L4191-04MS	15	1.13	9.8	9.4	95.4	
L4191-06	L4191-04MSD	16	1.13	9.8	9.4	95.4	
L4191-07	SS-31	17	1.13	9.82	8.88	89.2	
L4191-08	SS-30	18	1.15	9.55	8.93	92.6	
L4191-09	SS-29	19	1.1	9.84	9.36	94.5	
L4191-10	SS-28	20	1.13	9.56	8.61	88.7	
L4191-11	SS-27	21	1.16	9.97	9.66	96.5	



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 10/1/2020

OVENTEMP IN Celsius(°C): 107
Time IN: 17:00
In Date: 09/30/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:15
Out Date: 10/01/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB111222

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
L4191-12	SS-26	22	1.14	9.93	3.69	29.0	
L4191-13	SS-25	23	1.13	9.86	9.48	95.6	
L4191-14	SS-24	24	1.1	9.5	9.00	94.0	
L4191-15	SS-23	25	1.16	9.98	9.44	93.9	
L4191-16	SS-22	26	1.12	9.54	9.25	96.6	
L4191-17	SS-21	27	1.13	9.71	5.87	55.2	
L4191-18	SS-20	28	1.13	9.88	8.25	81.4	
L4191-19	SS-19	29	1.11	9.76	3.77	30.8	
L4191-20	SS-18	30	1.13	9.6	9.11	94.2	
L4191-21	SS-17	31	1.14	9.73	5.65	52.5	
L4191-22	SS-16	32	1.14	9.97	9.51	94.8	
L4192-01	RBS-BD001-01	33	1.00	2.00	2.00	100.0	Vegetative Debris,100% SOLIDS
L4193-01	B-12(3-4)	34	1.14	9.7	9.09	92.9	
L4193-02	B-12(9-11)	35	1.12	9.52	8.46	87.4	
L4193-03	B-12(13-14)	36	1.18	9.72	8.81	89.3	
L4193-04	B-13(11-13)	37	1.15	9.95	8.82	87.2	
L4193-05	B-14(9-11)	38	1.13	9.85	9.09	91.3	
L4193-06	B-14(11-13)	39	1.14	9.94	8.98	89.1	
L4193-07	B-14(13-14)	40	1.15	9.98	8.96	88.4	
L4200-01	PB-2020-WS-000-01	41	1.15	9.97	9.14	90.6	

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$