



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 10/2/2020

OVENTEMP IN Celsius(°C): 107
Time IN: 17:00
In Date: 10/01/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:16
Out Date: 10/02/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB111251

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
L3865-01	NB-08-093020	1	1.15	9.95	8.64	85.1	
L3865-02	NB-08-093020	2	1.15	9.96	8.66	85.2	
L3865-03	NB-08-093020-1	3	1.12	9.66	8.34	84.5	
L3865-05	NB-07-093020	4	1.15	9.95	8.76	86.5	
L3865-06	NB-07-093020	5	1.13	9.74	8.38	84.2	
L3865-07	NB-07-093020-1	6	1.18	9.75	8.4	84.2	
L3872-13	BU-03-093020	7	1.18	9.95	8.57	84.3	
L3872-14	BU-03-093020	8	1.16	9.98	8.93	88.1	
L3872-15	BU-03-093020-1	9	1.14	9.98	8.7	85.5	
L4195-01	RBR-200010	10	1.11	9.74	2.6	17.3	
L4195-02	RBR-200010	11	1.16	9.79	3.76	30.1	
L4195-03	RBR-251673	12	1.13	9.8	7.37	72.0	
L4195-04	RBR-251673	13	1.13	9.65	7.15	70.7	
L4196-01	RBR-251609	14	1.14	9.98	5.65	51.0	
L4196-02	RBR-251609	15	1.11	9.66	5.42	50.4	
L4213-01	B-15 (3-4)	16	1.13	9.93	9.06	90.1	
L4213-02	L4213-01MS	17	1.13	9.93	9.06	90.1	
L4213-03	L4213-01MSD	18	1.13	9.93	9.06	90.1	
L4213-04	B-15 (9-11)	19	1.14	9.88	9.16	91.8	
L4213-05	B-15 (13-14)	20	1.13	9.57	8.74	90.2	
L4213-06	B-16 (7-9)	21	1.13	9.82	9.11	91.8	
L4213-07	B-16 (9-11)	22	1.17	9.6	8.69	89.2	
L4213-08	B-16 (11-13)	23	1.13	9.7	8.5	86.0	
L4213-09	B-17 (3-4)	24	1.18	9.51	8.79	91.4	
L4213-10	B-17 (5-7)	25	1.14	9.64	9.19	94.7	
L4213-11	DUP-3	26	1.13	9.79	8.77	88.2	
L4213-12	B-17 (13-15)	27	1.13	9.83	8.81	88.3	
L4214-01	LOD-MDL-SOIL-01-QT4-2020	28	1.00	2.00	2.00	100.0	



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 10/2/2020

OVENTEMP IN Celsius(°C): 107
Time IN: 17:00
In Date: 10/01/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:16
Out Date: 10/02/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB111251

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
L4214-02	LOQ-SOIL-02-QT4-2020	29	1.00	2.00	2.00	100.0	
L4214-03	MDL-SOIL-03-QT4-2020	30	1.00	2.00	2.00	100.0	
L4214-04	LOD-MDL-SOIL-04-QT4-2020	31	1.00	2.00	2.00	100.0	
L4214-05	LOQ-SOIL-05-QT4-2020	32	1.00	2.00	2.00	100.0	
L4214-06	MDL-SOIL-06-QT4-2020	33	1.00	2.00	2.00	100.0	
L4219-01	CL-01-100120	63	1.14	9.83	9.06	91.1	
L4219-02	CL-01-100120	64	1.13	9.8	9.05	91.3	
L4219-03	CL-01-100120-1	65	1.15	9.98	9.1	90.0	
L4219-04	CL-01-100120-2	66	1.13	9.82	8.72	87.3	
L4219-05	CL-01-100120-3	67	1.19	9.78	8.81	88.7	
L4219-06	CL-01-100120-4	68	1.13	9.73	9.36	95.7	
L4219-07	CL-01-100120-5	69	1.15	9.96	8.59	84.4	
L4221-01	OK-02-100120	70	1.13	9.52	9.02	94.0	
L4221-02	OK-02-100120	71	1.17	9.68	8.7	88.5	
L4221-03	OK-02-100120-1	72	1.18	9.75	9.34	95.2	
L4221-04	OK-02-100120-2	73	1.17	9.76	9.54	97.4	
L4226-01	TR-05-100120	74	1.15	9.98	8.65	84.9	
L4226-02	TR-05-100120	75	1.1	9.74	8.49	85.5	
L4226-03	TR-05-100120-1	76	1.18	9.96	9.12	90.4	
L4228-01	PY-BIN-4	77	1.13	9.82	9.23	93.2	
L4228-02	PY-BIN-4	78	1.15	9.98	9.34	92.8	
L4230-01	RBR-29	79	1.13	9.65	9.38	96.8	
L4230-02	RBR-29	80	1.13	9.52	9.29	97.3	
L4230-03	26516	81	1.19	9.65	8.42	85.5	
L4230-04	26516	82	1.18	9.71	8.43	85.0	
L4230-05	2392	83	1.15	9.96	8.47	83.1	
L4230-06	2392	84	1.1	9.55	9.29	96.9	
L4230-07	24702	85	1.13	9.56	6.95	69.0	



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 10/2/2020

OVENTEMP IN Celsius(°C): 107
Time IN: 17:00
In Date: 10/01/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:16
Out Date: 10/02/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB111251

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
L4230-08	24702	86	1.16	9.96	7.69	74.2	
L4230-09	2015	87	1.13	9.82	8.44	84.1	
L4230-10	2015	88	1.15	9.97	9.14	90.6	
L4230-11	27080	89	1.00	2.00	2.00	100.0	OILY DEBRIS
L4230-12	27080	90	1.00	2.00	2.00	100.0	OILY DEBRIS
L4231-01	337	91	1.16	98.54	7.56	6.6	
L4231-02	337	92	1.1	9.85	7.98	78.6	
L4238-01	WASTE	34	1.13	9.79	8.42	84.2	
L4238-02	VOC	35	1.14	9.94	8.58	84.5	
L4238-03	1	36	1.18	9.73	8.39	84.3	
L4238-04	2	37	1.11	9.66	8.2	82.9	
L4238-05	3	38	1.1	9.8	8.44	84.4	
L4238-06	4	39	1.15	9.95	8.58	84.4	
L4238-07	5	40	1.14	9.98	8.64	84.8	
L4245-01	PE-27-(0.0-0.5)	41	1.14	9.96	8.52	83.7	
L4245-02	PE-27-(4.5-5.0)	42	1.18	9.66	8.2	82.8	
L4245-03	PE-27-(6.5-7.0)	43	1.15	9.95	9.16	91.0	
L4245-04	PE-28-(4.5-5.0)	44	1.18	9.55	8.58	88.4	
L4245-05	PE-28-(6.5-7.0)	45	1.1	9.81	8.78	88.2	
L4245-06	PE-29-(4.5-5.0)	46	1.18	9.97	9.21	91.4	
L4245-07	PE-29-(6.5-7.0)	47	1.12	9.5	8.31	85.8	
L4245-08	PE-30-(4.5-5.0)	48	1.14	9.96	9.8	98.2	
L4245-09	PE-30-(6.5-7.0)	49	1.15	9.92	9.16	91.3	
L4245-10	PE-31-(4.5-5.0)	50	1.14	9.9	9.00	89.7	
L4245-11	PE-31-(6.5-7.0)	51	1.16	9.77	8.8	88.7	
L4245-12	PE-32-(4.5-5.0)	52	1.1	9.82	8.89	89.3	
L4245-13	PE-32-(6.5-7.0)	53	1.13	9.97	9.84	98.5	
L4245-14	PE-33-(4.5-5.0)	54	1.14	9.98	8.7	85.5	



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 10/2/2020

OVENTEMP IN Celsius(°C): 107
Time IN: 17:00
In Date: 10/01/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:16
Out Date: 10/02/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB111251

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
L4245-15	PE-33- (6.5-7.0)	55	1.18	9.92	8.54	84.2	
L4245-16	PE-34- (4.5-5.0)	56	1.19	9.98	9.81	98.1	
L4245-17	PE-34- (6.5-7.0)	57	1.14	9.85	8.59	85.5	
L4245-18	PE-35- (6.5-7.0)	58	1.18	9.59	8.22	83.7	
L4245-19	PE-36- (6.5-7.0)	59	1.16	9.93	9.13	90.9	
L4245-20	PE-37- (6.5-7.0)	60	1.17	9.77	8.35	83.5	
L4245-21	PE-38- (6.5-7.0)	61	1.13	9.93	9.06	90.1	
L4245-22	DUP-02- (20201001)	62	1.16	9.58	8.34	85.3	
L4246-01	108885	93	1.00	2.00	2.00	100.0	OIL SAMPLE

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$