



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 10/5/2020

OVENTEMP IN Celsius(°C): 107
Time IN: 14:25
In Date: 10/02/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 07:50
Out Date: 10/03/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB111271

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
L4207-01	MBEQK5	1	1.15	9.69	3.62	28.9	
L4207-02	MBEQ27	2	1.15	9.98	4.34	36.1	
L4207-03	MBEQ28	3	1.14	9.64	3.27	25.1	
L4207-04	MBEQ29	4	1.15	9.66	3.23	24.4	
L4207-05	MBEQ30	5	1.14	9.96	6.11	56.3	
L4207-06	MBEQ31	6	1.18	9.72	3.81	30.8	
L4207-07	MBEQ32	7	1.14	9.94	3.44	26.1	
L4207-08	MBEQ32D	8	1.14	9.94	3.44	26.1	
L4207-09	MBEQ32S	9	1.14	9.94	3.44	26.1	
L4207-10	MBEQ33	10	1.11	9.97	5.17	45.8	
L4207-11	MBEQ34	11	1.14	9.98	5.83	53.1	
L4207-12	MBEQ35	12	1.14	9.92	3.72	29.4	
L4207-13	MBEQ36	13	1.18	9.98	4.33	35.8	
L4207-14	MBEQ71	14	1.12	9.61	6.97	68.9	
L4207-17	MBEQA7	15	1.15	9.98	7.02	66.5	
L4207-18	MBEQA8	16	1.15	9.95	7.91	76.8	
L4207-19	MBEQA9	17	1.15	9.93	7.92	77.1	

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$