

Analysis Method: 9034
 Parameter: Reactive Sulfide
 Run Number: LB111282

ANALYST: ketankumar
 SUPERVISOR REVIEW BY: apatel
 Constant: 16000
 Normality1: 0.025
 Normality2: 0.025

Reagent/Standard	Lot/Log #
Sodium thiosulfate titrant 0.025N	WP86633
IODINE SOLUTION .025N 1L	W2670
Starch Solution, 4L	W2703

Seq	Lab ID	True Value (mg/l)	DF	Initial Weight (g)	Final Volume (ml)	T1 (ml)	T2 Initial	T2 Final	T2 Diff. (ml)	T1 - T2 Diff (mL)	Value Corrected With Blank	Result (ppm)	Anal Date	Anal Time
1	PB132080BL		1	5.01	50	2.00	0.00	1.90	1.90	0.10	0.00	0.00	10/02/2020	15:10
2	L4190-04		1	5.02	50	2.00	1.90	3.74	1.84	0.16	0.06	4.78	10/02/2020	15:12
3	L4190-04DUP		1	5.04	50	2.00	3.74	5.60	1.86	0.14	0.04	3.17	10/02/2020	15:15
4	L4190-06		1	5.05	50	2.00	5.60	7.42	1.82	0.18	0.08	6.34	10/02/2020	15:18
5	L4190-12		1	5.04	50	2.00	0.00	1.88	1.88	0.12	0.02	1.59	10/02/2020	15:20
6	L4190-14		1	5.06	50	2.00	1.88	3.70	1.82	0.18	0.08	6.32	10/02/2020	15:22
7	L4190-18		1	5.02	50	2.00	3.70	5.50	1.80	0.20	0.10	7.97	10/02/2020	15:25
8	L4200-01		1	5.03	50	2.00	5.50	7.34	1.84	0.16	0.06	4.77	10/02/2020	15:28
9	L4219-08		1	5.04	50	2.00	0.00	1.82	1.82	0.18	0.08	6.35	10/02/2020	15:30
10	L4221-05		1	5.02	50	2.00	1.82	3.60	1.78	0.22	0.12	9.56	10/02/2020	15:32
11	L4226-04		1	5.03	50	2.00	3.60	5.40	1.80	0.20	0.10	7.95	10/02/2020	15:35

T1 = Titrant1

T2 = Titrant2

T2 Diff = T2 Final - T2 Initial

Value Corrected With Blank = ((T1 - T2 Diff) - Blank Correction(BL))

Result = ((T1 * Normality1) - ((T1 - Value Corrected With Blank) * Normality2)) * Constant / Initial Volume