



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 10/7/2020

OVENTEMP IN Celsius(°C): 107
Time IN: 17:10
In Date: 10/06/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:11
Out Date: 10/07/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB111328

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
L4221-06	OK-01-100620	58	1.15	9.75	9.33	95.1	
L4221-07	OK-01-100620	59	1.13	9.61	9.11	94.1	
L4221-08	OK-01-100620-1	60	1.13	9.76	9.36	95.4	
L4224-01	OR-02-100620	61	1.18	9.92	9.55	95.8	
L4224-02	OR-02-100620	62	1.15	9.75	9.4	95.9	
L4224-03	OR-02-100620-1	63	1.14	9.82	9.61	97.6	
L4224-04	OR-02-100620-2	64	1.19	9.78	9.22	93.5	
L4225-01	SU-04-100620	65	1.17	9.54	9.28	96.9	
L4225-02	SU-04-100620	66	1.19	9.76	9.58	97.9	
L4225-03	SU-04-100620-1	67	1.17	9.55	9.5	99.4	
L4272-01	MH-1- (6)	1	1.14	9.76	8.98	91.0	
L4272-02	MH-1- (6) -EPH-2	2	1.14	9.59	8.77	90.3	
L4272-03	MH-1- (6)	3	1.17	9.71	8.91	90.6	
L4272-04	MH-2- (6)	4	1.18	9.58	8.84	91.2	
L4272-05	MH-2- (6) -EPH-2	5	1.12	9.66	8.76	89.5	
L4272-06	MH-2- (6)	6	1.14	9.9	9.06	90.4	
L4272-07	MH-3- (6)	7	1.15	9.95	8.95	88.6	
L4272-08	MH-3- (6) -EPH-2	8	1.11	9.82	8.81	88.4	
L4272-09	MH-3- (6)	9	1.14	9.72	8.51	85.9	
L4272-10	GIS-1- (6)	10	1.18	9.68	8.6	87.3	
L4272-11	GIS-1- (6) -EPH-2	11	1.18	9.55	8.36	85.8	
L4272-12	GIS-1- (6)	12	1.15	9.97	8.33	81.4	
L4272-13	GIS-2- (6)	13	1.13	9.78	8.96	90.5	
L4272-14	GIS-2- (6) -EPH-2	14	1.14	9.51	8.7	90.3	
L4272-15	GIS-2- (6)	15	1.13	9.57	8.69	89.6	
L4272-16	GIS-3- (6)	16	1.14	9.62	9.16	94.6	
L4272-17	GIS-3- (6) -EPH-2	17	1.17	9.53	9.08	94.6	
L4272-18	GIS-3- (6)	18	1.11	9.79	9.32	94.6	



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 10/7/2020

OVENTEMP IN Celsius(°C): 107
Time IN: 17:10
In Date: 10/06/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:11
Out Date: 10/07/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB111328

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
L4272-19	GIS-4- (6)	19	1.16	9.75	9.26	94.3	
L4272-20	GIS-4- (6) -EPH-2	20	1.12	9.58	9.06	93.9	
L4272-21	GIS-4- (6)	21	1.14	9.83	9.33	94.2	
L4272-22	SC-1- (6.5)	22	1.13	9.58	8.7	89.6	
L4272-23	SC-1- (6.5) -EPH-2	23	1.16	9.71	8.78	89.1	
L4272-24	SC-1- (6.5)	24	1.13	9.59	8.65	88.9	
L4272-25	SC-2- (6.5)	25	1.14	9.98	9.16	90.7	
L4272-26	SC-2- (6.5) -EPH-2	26	1.1	9.86	9.02	90.4	
L4272-27	SC-2- (6.5)	27	1.18	9.66	8.86	90.6	
L4272-28	SC-3- (6.5)	28	1.18	9.75	8.85	89.5	
L4272-29	SC-3- (6.5) -EPH-2	29	1.13	9.88	9.05	90.5	
L4272-30	SC-3- (6.5)	30	1.14	9.54	8.7	90.0	
L4272-31	SC-4- (6.5)	31	1.14	9.94	9.26	92.3	
L4272-32	SC-4- (6.5) -EPH-2	32	1.17	9.68	8.83	90.0	
L4272-33	SC-4- (6.5)	33	1.18	9.56	8.7	89.7	
L4272-34	SC-5- (6.5)	34	1.14	9.66	8.79	89.8	
L4272-35	SC-5- (6.5) -EPH-2	35	1.13	9.83	8.98	90.2	
L4272-36	SC-5- (6.5)	36	1.14	9.97	9.04	89.5	
L4272-37	SC-6- (6.5)	37	1.18	9.69	9.14	93.5	
L4272-38	SC-6- (6.5) -EPH-2	38	1.18	9.66	9.09	93.3	
L4272-39	SC-6- (6.5)	39	1.13	9.77	9.15	92.8	
L4272-40	SC-7- (6.5)	40	1.17	9.97	9.24	91.7	
L4272-41	SC-7- (6.5) -EPH-2	41	1.16	9.75	8.99	91.2	
L4272-42	SC-7- (6.5)	42	1.14	9.65	8.84	90.5	
L4272-43	SC-8- (6.5)	43	1.19	9.77	9.33	94.9	
L4272-44	SC-8- (6.5) -EPH-2	44	1.16	9.82	9.4	95.2	
L4272-45	SC-8- (6.5)	45	1.14	9.89	9.48	95.3	
L4272-46	SC-9- (3.5)	46	1.18	9.9	9.08	90.6	



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 10/7/2020

OVENTEMP IN Celsius(°C): 107
Time IN: 17:10
In Date: 10/06/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:11
Out Date: 10/07/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB111328

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
L4272-47	SC-9-(3.5)-EPH-2	47	1.16	9.76	8.92	90.2	
L4272-48	SC-9-(3.5)	48	1.16	9.57	8.8	90.8	
L4272-49	SC-10-(3.5)	49	1.17	9.52	8.6	89.0	
L4272-50	SC-10-(3.5)-EPH-2	50	1.17	9.82	8.83	88.6	
L4272-51	SC-10-(3.5)	51	1.18	9.96	9.3	92.5	
L4272-52	SC-11-(3.5)	52	1.14	9.65	8.44	85.8	
L4272-53	SC-11-(3.5)-EPH-2	53	1.13	9.63	8.45	86.1	
L4272-54	SC-11-(3.5)	54	1.16	9.86	9.15	91.8	
L4272-55	SC-12-(3.5)	55	1.18	9.72	8.96	91.1	
L4272-56	SC-12-(3.5)-EPH-2	56	1.1	9.61	8.81	90.6	
L4272-57	SC-12-(3.5)	57	1.18	9.97	9.26	91.9	
L4279-01	S-1	68	1.12	9.74	4.35	37.5	
L4279-03	S-2	69	1.18	9.56	4.71	42.1	
L4281-01	23849	70	1.00	2.00	2.00	100.0	OILY DEBRIS
L4281-02	31193	71	1.13	9.59	8.89	91.7	
L4281-03	49-93	72	1.14	9.74	8.49	85.5	
L4282-01	31068	73	1.14	9.51	9.33	97.8	
L4282-02	989	74	1.00	2.00	2.00	100.0	OILY DEBRIS
L4282-03	68-89	75	1.13	9.53	9.36	98.0	

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$