



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 10/8/2020

OVENTEMP IN Celsius(°C): 107
Time IN: 17:20
In Date: 10/07/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:11
Out Date: 10/08/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB111366

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
L4291-01	PAV-B4-100620-0-1	1	1.18	9.69	8.84	90.0	
L4291-03	PAV-B4-100620-0-10	2	1.15	9.52	8.57	88.6	
L4291-05	PAV-B4-100620-10-20	3	1.12	9.81	9.00	90.7	
L4291-07	PAV-B4-100620-19-20	4	1.12	9.73	8.92	90.6	
L4291-09	PAV-B4-100620-DP	5	1.17	9.61	8.79	90.3	
L4291-11	PAV-B4-100620-DP	6	1.12	9.54	8.59	88.7	
L4292-01	PAV-B13-100620-0-1	7	1.14	9.9	9.22	92.2	
L4292-03	PAV-B13-100620-0-10	8	1.1	9.89	9.35	93.9	
L4292-05	PAV-B13-100620-15-16	9	1.15	9.96	8.98	88.9	
L4292-07	PAV-B13-100620-DP	10	1.18	9.58	8.99	93.0	
L4292-09	PAV-B13-100620-DP	11	1.12	9.83	9.33	94.3	
L4294-01	TP- (0-3) - (3-5) -COMP	30	1.12	9.55	9.04	94.0	
L4294-03	TP- (5) -GRAB	31	1.14	9.92	9.36	93.6	
L4294-05	TP- (0-3) - (3-5) -COMP-EPH-2	32	1.12	9.55	9.04	94.0	
L4295-01	DRUM-1	12	1.00	2.00	2.00	100.0	OILY DEBRIS
L4295-02	10-07A	13	1.00	2.00	2.00	100.0	PLASTIC
L4295-03	10-07B	14	1.00	2.00	2.00	100.0	PLASTIC
L4295-04	10-07C	15	1.00	2.00	2.00	100.0	PLASTIC
L4296-01	31281	16	1.00	2.00	2.00	100.0	OIL SAMPLE
L4297-01	RO-RBR-200031	33	1.13	9.62	3.39	26.6	Gel matrix
L4297-02	RO-RBR-200031	34	1.12	9.98	2.2	12.2	Gel matrix
L4297-03	RO-RBR-200034	35	1.15	9.98	1.98	9.4	Gel matrix
L4297-04	RO-RBR-200034	36	1.14	9.88	1.56	4.8	Gel matrix
L4298-01	RBR-559	17	1.12	9.71	8.88	90.3	
L4298-02	RBR-559	18	1.13	9.71	8.9	90.6	
L4300-01	231-232-COMP	19	1.14	9.55	6.81	67.4	
L4303-21	IDW-20201006	37	1.1	9.72	9.68	99.5	
L4306-01	OILY-DEBRIS-DRUM	38	1.00	2.00	2.00	100.0	OILY DEBRIS



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 10/8/2020

OVENTEMP IN Celsius(°C): 107
Time IN: 17:20
In Date: 10/07/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:11
Out Date: 10/08/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB111366

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
L4309-01	OWLS-HEAD-BH-02-A	20	1.18	9.72	9.21	94.0	
L4309-02	OWLS-HEAD-BH-02-B	21	1.14	9.59	9.28	96.3	
L4309-03	OWLS-HEAD-BH-04-A	22	1.18	9.64	9.13	94.0	
L4309-04	OWLS-HEAD-BH-04-B	23	1.14	9.6	9.1	94.1	
L4309-05	OWLS-HEAD-BH-05-A	24	1.18	9.97	9.26	91.9	
L4309-06	OWLS-HEAD-BH-05-B	25	1.13	9.86	8.9	89.0	
L4311-01	PAV-B10-100720-0-1	26	1.13	9.94	9.31	92.8	
L4311-03	PAV-B10-100720-0-10	27	1.1	9.69	8.91	90.9	
L4311-05	PAV-B10-100720-10-17	28	1.15	9.97	9.29	92.3	
L4311-07	PAV-B10-100720-16-17	29	1.1	9.71	9.16	93.6	

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$