



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 10/9/2020

OVENTEMP IN Celsius(°C): 107
Time IN: 14:30
In Date: 10/08/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:00
Out Date: 10/09/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB111373

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
L4267-01	BEQC3	1	1.14	9.94	4.48	38.0	
L4267-02	BEQC3MS	2	1.14	9.94	4.48	38.0	
L4267-03	BEQC3MSD	3	1.14	9.94	4.48	38.0	
L4267-04	BEQ55	4	1.18	9.97	6.17	56.8	
L4267-05	BEQ56	5	1.12	9.71	3.49	27.6	
L4267-06	BEQ57	6	1.12	9.55	4.56	40.8	
L4267-07	BEQ58	7	1.12	9.95	3.66	28.8	
L4267-08	BEQB9	8	1.12	9.66	5.54	51.8	
L4267-09	BEQC0	9	1.12	9.7	5.35	49.3	
L4267-10	BEQC1	10	1.18	9.66	5.04	45.5	
L4267-11	BEQC5	11	1.14	9.9	5.04	44.5	
L4267-12	BEQC6	12	1.14	9.92	5.65	51.4	
L4267-13	BEQC8	13	1.15	9.96	4.18	34.4	
L4267-14	BEQC9	14	1.13	9.66	5.21	47.8	
L4267-15	BEQ59	15	1.14	9.97	4.12	33.7	
L4267-16	BEQ60	16	1.18	9.93	3.39	25.3	
L4267-17	BEQ61	17	1.14	9.98	5.99	54.9	
L4267-18	BEQ62	18	1.13	9.62	4.08	34.7	
L4267-19	BEQ63	19	1.13	9.79	4.16	35.0	
L4267-20	BEQ64	20	1.16	9.98	5.55	49.8	
L4267-21	BEQ65	21	1.1	9.77	3.79	31.0	
L4267-22	BEQ66	22	1.15	9.97	5.86	53.4	

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$