



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 10/26/2020

OVENTEMP IN Celsius(°C): 107
Time IN: 16:30
In Date: 10/09/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:00
Out Date: 10/10/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB111407

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
L4225-05	SU-01-100920	22	1.18	9.83	9.05	91.0	
L4225-06	SU-01-100920	23	1.13	9.59	8.85	91.3	
L4225-07	SU-01-100920-1	24	1.18	9.64	8.88	91.0	
L4225-08	SU-01-100920-2	25	1.12	9.69	9.6	98.9	
L4335-01	01-A-01-B-01-C	1	1.00	2.00	2.00	100.0	CAULKING SAMPLE
L4336-05	SVOC-GPC-BLANK	2	1.00	2.00	2.00	100.0	
L4336-06	PEST-GPC-BLANK	3	1.00	2.00	2.00	100.0	
L4336-07	PEST-GPC-BLANK-SPIKE	4	1.00	2.00	2.00	100.0	
L4336-08	PCB-GPC-BLANK	5	1.00	2.00	2.00	100.0	
L4336-09	PCB-GPC-BLANK-SPIKE	6	1.00	2.00	2.00	100.0	
L4336-10	SVOC-GCP2-BLANK	7	1.00	2.00	2.00	100.0	
L4336-11	PEST-GCP2-BLANK	8	1.00	2.00	2.00	100.0	
L4336-12	PEST-GCP2-BLANK-SPIKE	9	1.00	2.00	2.00	100.0	
L4336-13	PCB-GCP2-BLANK	10	1.00	2.00	2.00	100.0	
L4336-14	PCB-GCP2-BLANK-SPIKE	11	1.00	2.00	2.00	100.0	
L4337-01	TP-4	17	1.15	9.96	8.66	85.2	
L4337-02	TP-4	18	1.14	9.98	8.79	86.5	
L4337-04	TP-4-EPH-2	19	1.17	9.98	8.58	84.1	
L4338-01	VNJ-235	29	1.14	9.78	9.51	96.9	
L4339-01	QUAD-1	30	1.00	2.00	2.00	100.0	CONCRETE SAMPLE
L4339-02	QUAD-2	31	1.00	2.00	2.00	100.0	CONCRETE SAMPLE
L4339-03	QUAD-3	32	1.00	2.00	2.00	100.0	CONCRETE SAMPLE
L4339-04	QUAD-4	33	1.00	2.00	2.00	100.0	CONCRETE SAMPLE
L4339-05	QUAD-5	34	1.00	2.00	2.00	100.0	CONCRETE SAMPLE
L4339-06	QUAD-6	35	1.00	2.00	2.00	100.0	CONCRETE SAMPLE
L4339-07	QUAD-7	36	1.00	2.00	2.00	100.0	CONCRETE SAMPLE
L4339-08	QUAD-8	37	1.00	2.00	2.00	100.0	CONCRETE SAMPLE
L4340-01	A-2332-MH11	20	1.13	9.68	9.29	95.4	



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 10/26/2020

OVENTEMP IN Celsius(°C): 107
Time IN: 16:30
In Date: 10/09/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:00
Out Date: 10/10/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB111407

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
L4340-03	A-2332-MH12	21	1.13	9.84	9.03	90.7	
L4341-01	STOCK-PILE	26	1.18	9.98	9.5	94.5	
L4341-02	STOCK-PILE	27	1.16	9.74	8.75	88.5	
L4341-02D	STOCK-PILEDUP	28	1.18	9.75	8.79	88.8	
L4342-06	BACKFILLOFFSITE-1-SO-1	12	1.18	9.6	8.61	88.2	
L4342-07	BACKFILLOFFSITE-2-SO-2	13	1.19	9.83	8.73	87.3	
L4342-08	BACKFILLOFFSITE-3-SO-3	14	1.13	9.58	8.66	89.1	
L4342-09	BACKFILLOFFSITE-3-SO-4	15	1.14	9.51	8.69	90.2	
L4342-10	BACKFILLOFFSITE-4-SO-5	16	1.13	9.91	8.86	88.0	

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$