



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 10/13/2020

OVENTEMP IN Celsius(°C): 106
Time IN: 15:00
In Date: 10/12/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:00
Out Date: 10/13/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB111434

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
L4285-01	BEQE6	1	1.15	9.67	3.52	27.8	
L4285-02	BEQE7	2	1.13	9.77	5.26	47.8	
L4285-03	BEQE8	3	1.15	9.98	7.63	73.4	
L4285-04	BEQE9	4	1.15	9.95	5.31	47.3	
L4285-05	BEQE9MS	5	1.15	9.95	5.31	47.3	
L4285-06	BEQE9MSD	6	1.15	9.95	5.31	47.3	
L4285-07	BEQF0	7	1.12	9.77	4.58	40.0	
L4285-08	BEQF1	8	1.13	9.74	4.57	40.0	
L4285-09	BEQG3	9	1.15	9.97	4.99	43.5	
L4285-11	BEQD7	10	1.13	9.82	4.62	40.2	
L4285-12	BEQD8	11	1.16	9.75	4.77	42.0	
L4285-13	BEQD9	12	1.18	9.98	7.48	71.6	
L4285-14	BEQE0	13	1.15	9.96	4.13	33.8	
L4285-15	BEQE1	14	1.13	9.66	4.88	44.0	
L4285-16	BEQE2	15	1.19	9.77	4.46	38.1	
L4285-17	BEQE3	16	1.14	9.97	7.04	66.8	
L4285-18	BEQE4	17	1.16	9.96	7.6	73.2	
L4285-19	BEQE5	18	1.13	9.97	7.9	76.6	

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$