



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 10/28/2020

OVENTEMP IN Celsius(°C): 106
Time IN: 15:00
In Date: 10/12/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:00
Out Date: 10/13/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB111435

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
L4320-01	MBEQE6	1	1.15	9.67	3.52	27.8	
L4320-02	MBEQE7	2	1.13	9.77	5.26	47.8	
L4320-03	MBEQE8	3	1.15	9.98	7.63	73.4	
L4320-04	MBEQE9	4	1.15	9.95	5.31	47.3	
L4320-05	MBEQE9D	5	1.15	9.95	5.31	47.3	
L4320-06	MBEQE9S	6	1.15	9.95	5.31	47.3	
L4320-07	MBEQF0	7	1.12	9.77	4.58	40.0	
L4320-08	MBEQF1	8	1.13	9.74	4.57	40.0	
L4320-09	MBEQG3	9	1.15	9.97	4.99	43.5	
L4320-11	MBEQD7	10	1.13	9.82	4.62	40.2	
L4320-12	MBEQD8	11	1.16	9.75	4.77	42.0	
L4320-13	MBEQD9	12	1.18	9.98	7.48	71.6	
L4320-14	MBEQE0	13	1.15	9.96	4.13	33.8	
L4320-15	MBEQE1	14	1.13	9.66	4.88	44.0	
L4320-16	MBEQE2	15	1.19	9.77	4.46	38.1	
L4320-17	MBEQE3	16	1.14	9.97	7.04	66.8	
L4320-18	MBEQE4	17	1.14	9.97	7.6	73.2	
L4320-19	MBEQE5	18	1.16	9.96	7.9	76.6	

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$