



PERCENT SOLID

Supervisor: apatel
Analyst: Alexander
Date: 10/15/2020

OVENTEMP IN Celsius(°C): 107
Time IN: 18:15
In Date: 10/14/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:10
Out Date: 10/15/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB111506

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
L4332-03	SHELL-L15-L14-TPH-2	1	1.19	9.78	6.73	64.5	
L4392-01	TP-1	2	1.16	9.61	7.8	78.6	
L4392-02	TP-1	3	1.16	9.61	7.8	78.6	
L4392-03	TP-1	4	1.16	9.61	7.8	78.6	
L4392-04	TP-1-EPH-2	5	1.14	9.53	7.63	77.4	
L4392-05	TP-2	6	1.16	9.59	8.66	89.0	
L4392-06	TP-2	7	1.16	9.59	8.66	89.0	
L4392-07	TP-2	8	1.16	9.59	8.66	89.0	
L4392-08	TP-2-EPH-2	9	1.15	9.54	8.7	90.0	
L4392-09	TP-3	10	1.15	9.63	8.99	92.5	
L4392-10	TP-3	11	1.15	9.63	8.99	92.5	
L4392-11	TP-3	12	1.15	9.63	8.99	92.5	
L4392-12	TP-3-EPH-2	13	1.14	9.58	7.56	76.1	
L4392-13	TP-4	14	1.15	9.84	8.53	84.9	
L4392-14	TP-4	15	1.15	9.84	8.53	84.9	
L4392-15	TP-4	16	1.15	9.84	8.53	84.9	
L4392-16	TP-4-EPH-2	17	1.13	9.93	8.77	86.8	
L4404-01	SHELL-L15-L14-NORTH	18	1.16	9.52	8.32	85.6	
L4404-02	SHELL-L15-L14-NORTH	19	1.18	9.67	8.36	84.6	
L4404-03	SHELL-L15-L14-NORTH-TPH-2	20	1.15	9.62	8.43	86.0	
L4406-01	JC-02-101420-A	21	1.15	9.67	8.74	89.1	
L4406-02	JC-02-101420-A	22	1.16	9.76	8.61	86.6	
L4406-03	JC-02-101420-1	23	1.16	9.7	9.21	94.3	
L4406-04	JC-02-101420-2	24	1.16	9.55	9.14	95.1	
L4406-05	JC-02-101420-3	25	1.16	9.56	8.94	92.6	
L4406-06	JC-02-101420-4	26	1.17	9.79	8.95	90.3	
L4406-07	JC-02-101420-5	27	1.16	9.53	8.51	87.8	
L4406-09	JC-02-101420-B	28	1.16	9.66	8.92	91.3	



PERCENT SOLID

Supervisor: apatel
Analyst: Alexander
Date: 10/15/2020

OVENTEMP IN Celsius(°C): 107
Time IN: 18:15
In Date: 10/14/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:10
Out Date: 10/15/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB111506

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
L4406-10	JC-02-101420-B	29	1.15	9.82	9.05	91.1	
L4406-11	JC-02-101420-1	30	1.15	9.74	9.19	93.6	
L4406-12	JC-02-101420-2	31	1.15	9.98	9.13	90.4	
L4407-01	FK-GF-02-062-A	32	1.14	9.52	8.7	90.2	
L4407-02	FK-GF-02-062-A	33	1.17	9.97	9.02	89.2	
L4407-03	FK-GF-02-062-B	34	1.18	9.56	8.58	88.3	
L4407-04	FK-GF-02-062-B	35	1.17	9.63	8.76	89.7	
L4408-02	T-1	36	1.18	9.78	8.11	80.6	
L4408-03	T-2	37	1.16	9.55	8.02	81.8	
L4408-04	T-3	38	1.13	9.72	7.86	78.3	
L4408-05	T-4	39	1.17	9.98	7.77	74.9	
L4408-06	T-5	40	1.16	9.96	9.1	90.2	
L4408-07	T-6	41	1.17	9.76	8.73	88.0	
L4408-08	T-7	42	1.15	9.55	8.00	81.5	
L4408-09	T-8	43	1.15	9.76	8.56	86.1	

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$