



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 10/16/2020

OVENTEMP IN Celsius(°C): 107
Time IN: 17:20
In Date: 10/15/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:11
Out Date: 10/16/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB111515

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
L4218-01	HD-01-101520	21	1.14	9.91	9.29	92.9	
L4218-02	HD-01-101520	22	1.14	9.72	8.47	85.4	
L4218-03	HD-01-101520-1	23	1.15	9.54	8.48	87.4	
L4219-09	CL-02-101520	24	1.14	9.63	8.95	92.0	
L4219-10	CL-02-101520	25	1.19	9.61	8.72	89.4	
L4219-11	CL-02-101520-1	26	1.12	9.51	9.00	93.9	
L4219-12	CL-02-101520-2	27	1.12	9.55	8.67	89.6	
L4219-13	CL-02-101520-3	28	1.15	9.96	9.26	92.1	
L4219-14	CL-02-101520-4	29	1.16	9.96	9.33	92.8	
L4219-15	CL-02-101520-5	30	1.1	9.82	9.03	90.9	
L4223-01	PL-01-101420	1	1.15	9.78	8.65	86.9	
L4223-02	PL-01-101420	2	1.18	9.88	8.51	84.3	
L4223-03	PL-01-101420-1	3	1.11	9.94	8.24	80.7	
L4223-04	PL-01-101420-2	4	1.13	9.81	9.32	94.4	
L4223-05	PL-01-101420-3	5	1.14	9.97	9.13	90.5	
L4223-06	PL-01-101420-4	6	1.16	9.97	8.77	86.4	
L4223-07	PL-01-101420-5	7	1.13	9.96	9.24	91.8	
L4224-06	OR-03-101520	31	1.17	9.94	9.42	94.1	
L4224-07	OR-03-101520	32	1.11	9.77	8.79	88.7	
L4224-08	OR-03-101520-1	33	1.12	9.55	8.36	85.9	
L4224-09	OR-03-101520-2	34	1.13	9.6	8.69	89.3	
L4227-01	BU-01-101420	8	1.13	9.86	8.83	88.2	
L4227-02	BU-01-101420	9	1.11	9.66	9.42	97.2	
L4227-03	BU-01-101420-1	10	1.14	9.98	9.19	91.1	
L4393-01	OILY-DEBRIS-DRUM	11	1.00	2.00	2.00	100.0	OILY DEBRIS
L4393-02	OILY-STONE-DRUMS	12	1.00	2.00	2.00	100.0	oily-stone
L4393-04	SP-COMP	13	1.12	9.88	8.17	80.5	
L4393-05	SP-TRUCK	14	1.17	9.9	8.47	83.6	



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 10/16/2020

OVENTEMP IN Celsius(°C): 107
Time IN: 17:20
In Date: 10/15/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:11
Out Date: 10/16/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB111515

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
L4394-01	BLUE-DRUM	15	1.00	2.00	2.00	100.0	OILY DEBRIS
L4395-01	OILY-DEBRIS-DRUMS-(1-4)	16	1.00	2.00	2.00	100.0	OILY DEBRIS
L4395-03	OILY-STONE-DRUMS-(1-2)	17	1.00	2.00	2.00	100.0	oily-stone
L4411-01	18618	35	1.00	2.00	2.00	100.0	OIL SAMPLE
L4411-02	32696	36	1.13	9.74	8.65	87.3	
L4413-01	30688	37	1.00	2.00	2.00	100.0	OILY DEBRIS
L4418-01	PB-2020-WS-000-02	18	1.14	9.94	8.72	86.1	
L4422-01	TP-3-(0-6FT)	19	1.13	9.75	8.39	84.2	
L4422-02	TP-(0-3)-(3-5)-COMP	20	1.12	9.55	9.04	94.0	
L4432-01	WCS-DD-G-01	38	1.14	9.97	8.52	83.6	

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$