



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 10/16/2020

OVENTEMP IN Celsius(°C): 107
Time IN: 15:10
In Date: 10/15/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:03
Out Date: 10/16/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB111526

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
L4402-01	C0AB4	1	1.15	9.94	7.69	74.4	
L4402-02	C0AB5	2	1.12	9.61	7.2	71.6	
L4402-03	C0AB6	3	1.13	9.72	7.58	75.1	
L4402-04	C0AB6MS	4	1.13	9.72	7.58	75.1	
L4402-05	C0AB6MSD	5	1.13	9.72	7.58	75.1	
L4402-06	C0AB7	6	1.13	9.72	7.53	74.5	
L4402-07	C0AB8	7	1.16	9.94	7.76	75.2	
L4402-08	C0AB9	8	1.18	9.73	5.98	56.1	
L4402-09	C0AC0	9	1.13	9.86	7.32	70.9	
L4402-10	C0AC1	10	1.15	9.94	6.93	65.8	
L4402-11	C0AC2	11	1.15	9.98	8.83	87.0	
L4402-12	C0AC3	12	1.11	9.94	6.6	62.2	
L4402-13	C0AC4	13	1.15	9.92	8.83	87.6	
L4402-14	C0AC5	14	1.13	9.9	6.44	60.5	
L4402-15	C0AC7	15	1.15	9.95	5.72	51.9	
L4402-16	C0AC8	16	1.1	9.81	5.16	46.6	
L4402-17	C0AC9	17	1.16	9.94	5.96	54.7	
L4402-18	C0AD0	18	1.13	9.64	5.58	52.3	

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$