



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 10/23/2020

OVENTEMP IN Celsius(°C): 107
Time IN: 17:00
In Date: 10/22/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:02
Out Date: 10/23/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB111651

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
L4485-01	MDL-SOIL-01	1	1.00	2.00	2.00	100.0	
L4485-02	MDL-SOIL-02	2	1.00	2.00	2.00	100.0	
L4485-03	MDL-SOIL-03	3	1.00	2.00	2.00	100.0	
L4485-04	MDL-SOIL-04	4	1.00	2.00	2.00	100.0	
L4485-05	MDL-SOIL-05	5	1.00	2.00	2.00	100.0	
L4485-06	MDL-SOIL-06	6	1.00	2.00	2.00	100.0	
L4485-07	MDL-SOIL-07	7	1.00	2.00	2.00	100.0	
L4485-15	MDL-ME-SOIL-01	8	1.00	2.00	2.00	100.0	
L4485-16	MDL-ME-SOIL-02	9	1.00	2.00	2.00	100.0	
L4485-17	MDL-ME-SOIL-03	10	1.00	2.00	2.00	100.0	
L4485-18	MDL-ME-SOIL-04	11	1.00	2.00	2.00	100.0	
L4485-19	MDL-ME-SOIL-05	12	1.00	2.00	2.00	100.0	
L4485-20	MDL-ME-SOIL-06	13	1.00	2.00	2.00	100.0	
L4485-21	MDL-ME-SOIL-07	14	1.00	2.00	2.00	100.0	
L4486-01	MS-SOIL	26	1.12	9.88	9.09	91.0	
L4486-02	MS-SOIL	27	1.12	9.85	8.94	89.6	
L4487-01	D-85	28	1.00	2.00	2.00	100.0	OIL SAMPLE
L4487-02	D-88	29	1.00	2.00	2.00	100.0	OIL SAMPLE
L4487-03	D-89	30	1.00	2.00	2.00	100.0	OIL SAMPLE
L4487-04	D-90	31	1.00	2.00	2.00	100.0	OIL SAMPLE
L4487-05	D-91	32	1.00	2.00	2.00	100.0	OIL SAMPLE
L4487-06	D-92	33	1.00	2.00	2.00	100.0	OIL SAMPLE
L4488-01	1618	34	1.00	2.00	2.00	100.0	OIL SAMPLE
L4488-04	202	35	1.15	9.79	9.36	95.0	
L4488-05	202	36	1.15	9.96	9.05	89.7	
L4491-01	COMP-1	15	1.14	9.67	7.94	79.7	
L4491-02	COMP-2	16	1.13	9.86	7.9	77.5	
L4491-03	COMP-3	17	1.16	9.67	8.87	90.6	



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 10/23/2020

OVENTEMP IN Celsius(°C): 107
Time IN: 17:00
In Date: 10/22/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:02
Out Date: 10/23/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB111651

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
L4494-01	AB-EH	19	1.18	9.96	9.16	90.9	
L4494-02	AB-EN	20	1.12	9.64	9.29	95.9	
L4494-03	AB-EH	21	1.18	9.96	9.16	90.9	
L4494-04	AB-EN	22	1.12	9.64	9.29	95.9	
L4495-01	63202	37	1.00	2.00	2.00	100.0	OIL SAMPLE
L4496-01	SA-01-102220	23	1.18	9.62	8.55	87.3	
L4496-02	SA-01-102220	24	1.16	9.86	8.9	89.0	
L4496-03	SA-01-102220-1	25	1.11	9.72	8.75	88.7	
L4498-01	38092-121119	18	1.00	2.00	2.00	100.0	P.E.SAMPLE

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$