



PERCENT SOLID

Supervisor: apatel
Analyst: ketankumar
Date: 11/5/2020

OVENTEMP IN Celsius(°C): 107
Time IN: 17:00
In Date: 10/26/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 108
Time OUT: 09:03
Out Date: 10/27/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB111688

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
L4512-01	WASTE	1	1.12	9.78	8.8	88.7	
L4512-02	VOC	2	1.14	9.56	8.66	89.3	
L4512-03	1	3	1.18	9.61	8.57	87.7	
L4512-04	2	4	1.13	9.63	8.61	88.0	
L4512-05	3	5	1.15	9.94	9.11	90.6	
L4512-06	4	6	1.16	9.97	9.42	93.8	
L4512-07	5	7	1.13	9.67	8.65	88.1	
L4513-01	WASTE	8	1.13	9.73	8.76	88.7	
L4513-02	VOC	9	1.14	9.86	9.03	90.5	
L4513-03	1	10	1.14	9.93	9.18	91.5	
L4513-04	2	11	1.17	9.77	8.78	88.5	
L4513-05	3	12	1.15	9.97	9.33	92.7	
L4513-06	4	13	1.12	9.74	8.75	88.5	
L4513-07	5	14	1.13	9.83	8.85	88.7	
L4514-01	102220HU01-P	15	1.00	2.00	2.00	100.0	CAULKING SAMPLE
L4515-01	OS-1-BLDG-9-N-ELEVATOR	16	1.00	2.00	2.00	100.0	OIL SAMPLE
L4515-02	OS-2-BLDG-9-S-ELEVATOR	17	1.00	2.00	2.00	100.0	OIL SAMPLE
L4515-03	OS-3-BLDG-7-ELEVATOR	18	1.00	2.00	2.00	100.0	OIL SAMPLE
L4515-04	OS-4-BLDG-6-ELEVATOR	19	1.00	2.00	2.00	100.0	OIL SAMPLE
L4515-05	OS-5-BLDG-5-ELEVATOR	20	1.00	2.00	2.00	100.0	OIL SAMPLE
L4515-06	OS-6-BLDG-4-ELEVATOR	21	1.00	2.00	2.00	100.0	OIL SAMPLE
L4515-07	OS-7-BLDG-3-ELEVATOR	22	1.00	2.00	2.00	100.0	OIL SAMPLE
L4517-01	EFB666Y-END-1	23	1.00	2.00	2.00	100.0	CABLE SAMPLE
L4517-02	EFB666Y-END-2	24	1.00	2.00	2.00	100.0	CABLE SAMPLE
L4517-03	WOOD-END-1	25	1.00	2.00	2.00	100.0	CABLE SAMPLE
L4517-04	WOOD-END-2	26	1.00	2.00	2.00	100.0	CABLE SAMPLE
L4517-05	KMA858E-END-1	27	1.00	2.00	2.00	100.0	CABLE SAMPLE
L4517-06	KMA858E-END-2	28	1.00	2.00	2.00	100.0	CABLE SAMPLE



PERCENT SOLID

Supervisor: apatel
Analyst: ketankumar
Date: 11/5/2020

OVENTEMP IN Celsius(°C): 107
Time IN: 17:00
In Date: 10/26/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 108
Time OUT: 09:03
Out Date: 10/27/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB111688

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
L4517-07	BC137121LOK-END-1	29	1.00	2.00	2.00	100.0	CABLE SAMPLE
L4517-08	BC137121LOK-END-2	30	1.00	2.00	2.00	100.0	CABLE SAMPLE
L4518-01	HEE0841-END-1	31	1.00	2.00	2.00	100.0	CABLE SAMPLE
L4518-02	HEE0841-END-2	32	1.00	2.00	2.00	100.0	CABLE SAMPLE
L4519-01	BC270885LOK-END-1	33	1.00	2.00	2.00	100.0	CABLE SAMPLE
L4519-02	BC270885LOK-END-1	34	1.00	2.00	2.00	100.0	CABLE SAMPLE
L4521-01	OIL	41	1.00	2.00	2.00	100.0	OIL SAMPLE
L4521-05	S-1	42	1.12	9.7	8.27	83.3	
L4521-06	S-2	43	1.14	9.97	7.9	76.6	
L4521-07	S-3	44	1.16	9.71	8.46	85.4	
L4521-08	S-4	45	1.14	9.94	7.87	76.5	
L4521-09	S-5	46	1.12	9.55	7.99	81.5	
L4521-10	DUP	47	1.17	9.96	8.22	80.2	
L4523-01	CS-37	35	1.13	9.78	9.42	95.8	
L4523-02	CS-38	36	1.11	9.79	9.46	96.2	
L4523-03	CS-39	37	1.18	9.62	9.29	96.1	
L4523-04	CS-40	38	1.18	9.97	9.61	95.9	
L4523-05	CS-41	39	1.15	9.96	9.66	96.6	
L4523-06	CS-42	40	1.19	9.93	9.61	96.3	

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$