

Analysis Method: 9034
 Parameter: Reactive Sulfide
 Run Number: LB111719

ANALYST: ketankumar
 SUPERVISOR REVIEW BY: apatel
 Constant: 16000
 Normality1: 0.025
 Normality2: 0.025

Reagent/Standard	Lot/Log #
Sodium thiosulfate titrant 0.025N	WP87313
IODINE SOLUTION .025N 1L	W2670
Starch Solution, 4L	W2703

Seq	Lab ID	True Value (mg/l)	DF	Initial Weight (g)	Final Volume (ml)	T1 (ml)	T2 Initial	T2 Final	T2 Diff. (ml)	T1 - T2 Diff (mL)	Value Corrected With Blank	Result (ppm)	Anal Date	Anal Time
1	PB132655BL		1	5.02	50	2.00	0.00	1.90	1.90	0.10	0.00	0.00	10/28/2020	11:32
2	L4219-24		1	5.02	50	2.00	1.90	3.76	1.86	0.14	0.04	3.19	10/28/2020	11:18
3	L4219-24DUP		1	5.04	50	2.00	3.76	5.64	1.88	0.12	0.02	1.59	10/28/2020	11:22
4	L4224-15		1	5.05	50	2.00	5.64	7.48	1.84	0.16	0.06	4.75	10/28/2020	11:25
5	L4526-02		1	5.05	50	2.00	0.00	1.82	1.82	0.18	0.08	6.34	10/28/2020	11:28
6	L4527-02		1	5.03	50	2.00	1.82	3.66	1.84	0.16	0.06	4.77	10/28/2020	11:32
7	L4527-04		1	5.06	50	2.00	3.66	5.52	1.86	0.14	0.04	3.16	10/28/2020	11:35
8	L4527-06		1	5.03	50	2.00	5.52	7.32	1.80	0.20	0.10	7.95	10/28/2020	11:38
9	L4528-02		1	5.02	50	2.00	0.00	1.84	1.84	0.16	0.06	4.78	10/28/2020	11:40
10	L4528-04		1	5.00	50	2.00	1.84	3.70	1.86	0.14	0.04	3.20	10/28/2020	11:42

T1 = Titrant1

T2 = Titrant2

T2 Diff = T2 Final - T2 Initial

Value Corrected With Blank = ((T1 - T2 Diff) - Blank Correction(BL))

Result = ((T1 * Normality1) - ((T1 - Value Corrected With Blank) * Normality2)) * Constant / Initial Volume