



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 10/30/2020

OVENTEMP IN Celsius(°C): 107
Time IN: 16:00
In Date: 10/29/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:00
Out Date: 10/30/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB111751

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
L4556-13	MX2Q03	1	1.00	2.00	2.00	100.0	p.e. sample
L4556-14	MX2Q03D	2	1.00	2.00	2.00	100.0	p.e. sample
L4556-15	MX2Q03S	3	1.00	2.00	2.00	100.0	p.e. sample
L4556-16	MX2Q04	4	1.00	2.00	2.00	100.0	p.e. sample
L4556-17	MX2Q04D	5	1.00	2.00	2.00	100.0	p.e. sample
L4556-18	MX2Q04S	6	1.00	2.00	2.00	100.0	p.e. sample
L4556-19	MX2Q05	7	1.00	2.00	2.00	100.0	p.e. sample
L4556-20	MX2Q05D	8	1.00	2.00	2.00	100.0	p.e. sample
L4556-21	MX2Q05S	9	1.00	2.00	2.00	100.0	p.e. sample
L4556-22	MX2Q06	10	1.00	2.00	2.00	100.0	p.e. sample
L4556-23	MX2Q06D	11	1.00	2.00	2.00	100.0	p.e. sample
L4556-24	MX2Q06S	12	1.00	2.00	2.00	100.0	p.e. sample

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$