



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 11/2/2020

OVENTEMP IN Celsius(°C): 107
Time IN: 16:25
In Date: 10/30/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:11
Out Date: 10/31/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB111775

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
L4570-01	MW-PIER-S	1	1.00	2.00	2.00	100.0	CONCRETE SAMPLE
L4570-02	MW-FOUND-T1	2	1.00	2.00	2.00	100.0	CONCRETE SAMPLE
L4572-05	SVOC-GPC-BLANK	3	1.00	2.00	2.00	100.0	
L4572-06	PEST-GPC-BLANK	4	1.00	2.00	2.00	100.0	
L4572-07	PEST-GPC-BLANK-SPIKE	5	1.00	2.00	2.00	100.0	
L4572-08	PCB-GPC-BLANK	6	1.00	2.00	2.00	100.0	
L4572-09	PCB-GPC-BLANK-SPIKE	7	1.00	2.00	2.00	100.0	
L4572-10	SVOC-GCP2-BLANK	8	1.00	2.00	2.00	100.0	
L4572-11	PEST-GCP2-BLANK	9	1.00	2.00	2.00	100.0	
L4572-12	PEST-GCP2-BLANK-SPIKE	10	1.00	2.00	2.00	100.0	
L4572-13	PCB-GCP2-BLANK	11	1.00	2.00	2.00	100.0	
L4572-14	PCB-GCP2-BLANK-SPIKE	12	1.00	2.00	2.00	100.0	
L4579-01	SS-35	13	1.14	9.66	9.03	92.6	
L4579-02	SS-36	14	1.13	9.76	8.71	87.8	
L4579-03	SS-37	15	1.13	9.82	8.85	88.8	
L4579-04	SS-38	16	1.14	9.98	8.88	87.6	
L4579-05	SS-39	17	1.16	9.77	8.6	86.4	
L4579-06	SS-40	18	1.13	9.71	8.9	90.6	
L4579-07	SS-41	19	1.12	9.88	9.14	91.6	
L4579-08	SS-42	20	1.15	9.83	8.9	89.3	
L4579-09	SS-43	21	1.13	9.67	8.72	88.9	
L4579-10	SS-44	22	1.16	9.97	9.22	91.5	
L4579-11	SS-45	23	1.15	9.8	8.81	88.6	
L4579-12	SS-46	24	1.15	9.96	8.93	88.3	
L4579-13	SS-47	25	1.15	9.95	9.00	89.2	
L4579-14	L4579-13MS	26	1.15	9.95	9.00	89.2	
L4579-15	L4579-13MSD	27	1.15	9.95	9.00	89.2	
L4579-16	SS-48	28	1.16	9.98	9.37	93.1	



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 11/2/2020

OVENTEMP IN Celsius(°C): 107
Time IN: 16:25
In Date: 10/30/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:11
Out Date: 10/31/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB111775

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
L4579-17	SS-98	29	1.14	9.68	8.84	90.2	
L4579-18	SS-49	30	1.16	9.97	9.39	93.4	
L4579-19	SS-50	31	1.13	9.58	8.62	88.6	
L4579-20	SS-51	32	1.16	9.74	9.08	92.3	
L4579-21	SS-52	33	1.14	9.77	8.94	90.4	
L4579-22	SS-53	34	1.12	9.52	8.18	84.0	
L4580-01	SS-54	35	1.13	9.85	8.51	84.6	
L4580-02	SS-104	36	1.16	9.95	9.31	92.7	
L4580-03	SS-55	37	1.12	9.74	8.67	87.6	
L4580-04	L4580-03MS	38	1.12	9.74	8.67	87.6	
L4580-05	L4580-03MSD	39	1.12	9.74	8.67	87.6	
L4580-06	SS-56	40	1.14	9.96	8.76	86.4	
L4580-07	SS-57	41	1.14	9.8	8.9	89.6	
L4580-08	SS-58	42	1.14	9.93	8.77	86.8	
L4580-09	SS-59	43	1.13	9.91	8.53	84.3	
L4580-10	SS-60	44	1.14	9.56	8.52	87.6	
L4580-11	SS-61	45	1.15	9.94	8.36	82.0	
L4580-12	SS-62	46	1.16	9.9	8.17	80.2	
L4580-13	SS-63	47	1.15	9.87	8.53	84.6	
L4580-14	SS-64	48	1.14	9.95	8.99	89.1	
L4580-15	SS-65	49	1.17	9.93	9.22	91.9	
L4580-16	SS-66	50	1.16	9.69	8.52	86.3	
L4580-17	SS-67	51	1.15	9.9	9.42	94.5	
L4580-18	SS-68	52	1.15	9.91	8.82	87.6	

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$