



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 11/3/2020

OVENTEMP IN Celsius(°C): 107
Time IN: 17:15
In Date: 11/02/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:22
Out Date: 11/03/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB111793

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
L4600-01	CRT-17	9	1.14	9.69	8.57	86.9	
L4600-05	CRT-23	10	1.14	9.94	8.91	88.3	
L4600-09	CRT-19	11	1.15	9.95	9.36	93.3	
L4600-13	CRT-6	12	1.16	9.74	8.69	87.8	
L4600-20	CRT-17-EPH-2	13	1.11	9.51	8.54	88.5	
L4600-21	CRT-23-EPH-2	14	1.11	9.74	8.77	88.8	
L4600-22	CRT-19-EPH-2	15	1.17	9.98	9.39	93.3	
L4600-23	CRT-6-EPH-2	16	1.16	9.83	8.84	88.6	
L4605-01	CRT-01	17	1.16	9.94	9.16	91.1	
L4605-09	CRT-01-EPH-2	18	1.14	9.94	9.4	93.9	
L4606-01	S-1	1	1.15	9.6	8.66	88.9	
L4606-02	S-1A	2	1.16	9.77	8.9	89.9	
L4608-01	OIL-1	3	1.00	2.00	2.00	100.0	OIL SAMPLE
L4608-02	PE-1	4	1.19	9.89	7.41	71.5	
L4608-03	PE-2	5	1.15	9.95	6.88	65.1	
L4608-04	PE-3	6	1.13	9.76	7.45	73.2	
L4608-05	PE-4	7	1.13	9.93	7.16	68.5	
L4608-06	DUP	8	1.17	9.98	7.4	70.7	
L4609-01	12016	19	1.14	9.98	8.6	84.4	
L4609-02	12016-EPH-2	20	1.15	9.94	8.93	88.5	

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$