



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 11/4/2020

OVENTEMP IN Celsius(°C): 107
Time IN: 17:15
In Date: 11/03/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:15
Out Date: 11/04/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB111811

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
L4591-01	OK-01-110320	27	1.16	9.83	9.26	93.4	
L4591-02	OK-01-110320	28	1.12	9.58	9.04	93.6	
L4596-01	TR-05-110320	31	1.18	9.72	9.03	91.9	
L4596-02	TR-05-110320	32	1.18	9.97	9.14	90.6	
L4611-01	BW-SOIL	1	1.14	9.61	8.57	87.7	
L4613-01	VNJ-219	33	1.14	9.94	9.13	90.8	
L4613-02	VNJ-219	34	1.14	9.94	9.13	90.8	
L4613-03	282	35	1.19	9.79	9.14	92.4	
L4613-04	282	36	1.19	9.79	9.14	92.4	
L4617-01	PE-1- (0.5-1.0)	2	1.15	9.84	8.43	83.8	
L4617-02	PE-2- (0.5-1.0)	3	1.14	9.73	7.95	79.3	
L4617-03	PE-3- (0.5-1.0)	4	1.13	9.69	7.86	78.6	
L4617-04	PE-4- (0.5-1.0)	5	1.15	9.67	7.58	75.5	
L4617-05	PE-5- (0.5-1.0)	6	1.13	9.97	8.52	83.6	
L4617-06	DUPLICATE	7	1.16	9.96	8.58	84.3	
L4617-10	T-1	8	1.00	2.00	2.00	100.0	OIL SAMPLE
L4619-01	TP-2-COMP	9	1.18	9.9	8.73	86.6	
L4619-02	TP-2-EPH-2.5-3.0	10	1.13	9.74	8.72	88.2	
L4619-03	TP2-EPH-6.0-6.5	11	1.15	9.96	8.68	85.5	
L4619-04	TP-1-3-4-COMP	12	1.13	9.97	9.26	92.0	
L4619-05	TP-1-3-4-EPH-2.5-3.0	13	1.11	9.61	8.58	87.9	
L4619-06	TP-1-3-4-EPH-4.5-5.0	14	1.14	9.97	8.59	84.4	
L4621-01	PAV-B8-110220-0-1	15	1.14	9.77	8.74	88.1	
L4621-03	PAV-B8-110320-0-10	16	1.1	9.69	8.74	88.9	
L4621-05	PAV-B8-110320-10-20	17	1.14	9.65	8.29	84.0	
L4621-07	PAV-B8-110320-19-20	18	1.13	9.75	8.7	87.8	
L4622-09	PAV-B7-110220-0-1	29	1.13	9.76	8.69	87.6	
L4622-11	PAV-B7-110220-DP	30	1.1	9.63	8.57	87.6	



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 11/4/2020

OVENTEMP IN Celsius(°C): 107
Time IN: 17:15
In Date: 11/03/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:15
Out Date: 11/04/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB111811

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
L4623-01	WASTE	19	1.15	9.96	8.76	86.4	
L4623-02	VOC	20	1.15	9.61	8.35	85.1	
L4623-03	1	21	1.14	9.95	9.04	89.7	
L4623-04	2	22	1.11	9.56	8.55	88.0	
L4623-05	3	23	1.14	9.95	8.91	88.2	
L4623-06	4	24	1.12	9.83	8.82	88.4	
L4623-07	5	25	1.14	9.77	8.65	87.0	
L4623-08	6	26	1.13	9.82	8.82	88.5	

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$