



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 11/12/2020

OVENTEMP IN Celsius(°C): 107
Time IN: 17:00
In Date: 11/11/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 07:42
Out Date: 11/12/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB111920

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
L4587-01	HR-02-111120	11	1.13	9.75	9.29	94.7	
L4587-02	HR-02-111120-EPH-2	12	1.14	9.67	9.25	95.1	
L4595-01	SU-03-111120	19	1.13	9.57	9.1	94.4	
L4595-02	SU-03-111120	20	1.16	9.71	9.12	93.1	
L4683-01	OILY-DEBRIS-01-02-03	1	1.00	2.00	2.00	100.0	OILY-DEBRIS
L4683-05	TOTE	2	1.00	2.00	2.00	100.0	OIL SAMPLE
L4683-06	OILY-WATER-001-002-003	3	1.00	2.00	2.00	100.0	OIL SAMPLE
L4685-01	72-11978	4	1.12	9.72	8.44	85.1	
L4685-02	72-11978	5	1.16	9.92	8.38	82.4	
L4699-01	STOCKPILE	6	1.13	9.77	4.31	36.8	
L4702-01	M-5 (0-10)	7	1.13	9.8	9.08	91.7	
L4702-02	M-5 (10-20)	8	1.16	9.96	8.76	86.4	
L4702-03	M-6 (0-10)	9	1.1	9.65	8.43	85.7	
L4702-04	M-99 (0-10)	10	1.13	9.94	8.59	84.7	
L4704-01	KL-SOIL	35	1.15	9.94	8.55	84.2	
L4704-02	KL-SOIL-EPH-2	36	1.13	9.67	7.99	80.3	
L4705-01	COMP-1	21	1.19	9.77	9.52	97.1	
L4705-03	COMP-2	22	1.16	9.71	9.6	98.7	
L4705-05	32785	23	1.00	2.00	2.00	100.0	plastic sample
L4706-01	COMP-30706-30707-40854	24	1.00	2.00	2.00	100.0	oily rags
L4706-02	16874	25	1.1	9.88	9.7	97.9	
L4714-01	EROC-SP-1	13	1.14	9.55	8.02	81.8	
L4714-02	EROC-SP-2	14	1.13	9.56	7.92	80.5	
L4715-01	M-6 (10-20)	37	1.15	9.97	9.06	89.7	
L4716-01	STOCKPILE-1	38	1.16	9.94	8.73	86.2	
L4717-01	R-LA-35-3-WC	15	1.13	9.58	8.21	83.8	
L4717-02	R-LA-35-3	16	1.13	9.52	8.17	83.9	
L4718-01	R-LA-38-1-WC	17	1.15	9.98	8.24	80.3	



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 11/12/2020

OVENTEMP IN Celsius(°C): 107
Time IN: 17:00
In Date: 11/11/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 07:42
Out Date: 11/12/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB111920

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
L4718-02	R-LA-38-1	18	1.15	9.94	8.09	79.0	
L4719-01	TP-3-S1	26	1.13	9.65	8.61	87.8	
L4719-02	TP-3-S2	27	1.16	9.53	8.46	87.2	
L4719-03	TP-3-S	28	1.17	9.64	8.65	88.3	
L4719-05	TP-3-D1	29	1.1	9.54	8.98	93.4	
L4719-06	TP-3-D2	30	1.16	9.81	9.14	92.3	
L4719-07	TP-3-D	31	1.17	9.76	9.16	93.0	
L4719-09	TP-3A-S1	32	1.19	9.66	8.53	86.7	
L4719-10	TP-3A-S2	33	1.15	9.68	8.8	89.7	
L4719-11	TP-3A-S	34	1.11	9.79	8.87	89.4	

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$