



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 11/12/2020

OVENTEMP IN Celsius(°C): 107
Time IN: 14:00
In Date: 11/11/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 07:45
Out Date: 11/12/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB111923

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
L4708-01	MA42S3	1	1.15	9.58	8.87	91.6	
L4708-02	MA42S4	2	1.13	9.95	8.61	84.8	
L4708-03	MA42S5	3	1.12	9.8	8.66	86.9	
L4708-04	MA42S6	4	1.13	9.7	7.92	79.2	
L4708-05	MA42S7	5	1.14	9.62	9.11	94.0	
L4708-06	MA42S8	6	1.11	9.96	9.55	95.4	
L4708-07	MA42S9	7	1.13	9.86	9.47	95.5	
L4708-08	MA42S9D	8	1.13	9.86	9.47	95.5	
L4708-09	MA42S9S	9	1.13	9.86	9.47	95.5	
L4708-10	MA42T0	10	1.15	9.63	9.08	93.5	
L4708-11	MA42T1	11	1.14	9.54	8.93	92.7	
L4708-12	MA42T2	12	1.13	9.54	9.22	96.2	
L4708-13	MA42T3	13	1.14	9.9	8.13	79.8	
L4708-14	MA42T4	14	1.12	9.9	8.8	87.5	
L4708-15	MA42T5	15	1.15	9.93	8.09	79.0	
L4708-16	MA42T6	16	1.13	9.98	8.45	82.7	
L4708-17	MA42T7	17	1.14	9.7	8.55	86.6	
L4708-18	MA42T8	18	1.11	9.72	7.94	79.3	
L4708-19	MA42T9	19	1.14	9.68	8.29	83.7	
L4708-20	MA42W0	20	1.13	9.62	7.48	74.8	
L4708-21	MA42W1	21	1.00	2.00	2.00	100.0	P.E. SAMPLE

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$