



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 11/30/2020

OVENTEMP IN Celsius(°C): 107
Time IN: 14:10
In Date: 11/13/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 07:42
Out Date: 11/14/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB111965

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
L4709-01	MA42W5	1	1.15	9.62	8.16	82.8	
L4709-02	MA42W6	2	1.14	9.95	8.9	88.1	
L4709-03	MA42W7	3	1.12	9.66	8.9	91.1	
L4709-04	MA42W7D	4	1.13	9.66	8.9	91.1	
L4709-05	MA42W7S	5	1.13	9.66	8.9	91.1	
L4709-06	MA42W8	6	1.13	9.76	8.53	85.7	
L4709-07	MA42W9	7	1.14	9.67	7.97	80.1	
L4709-08	MA42X1	8	1.15	9.9	9.09	90.7	
L4709-09	MA42X2	9	1.14	9.91	8.49	83.8	
L4709-10	MA42X3	10	1.15	9.6	8.45	86.4	
L4709-11	MA42X5	11	1.14	9.55	8.58	88.5	
L4709-12	MA42X6	12	1.15	9.74	8.69	87.8	
L4709-13	MA42X7	13	1.14	9.88	9.16	91.8	
L4709-14	MA42X8	14	1.12	9.86	8.84	88.3	
L4709-15	MA42X9	15	1.14	9.98	9.19	91.1	
L4709-16	MA42Y1	16	1.14	9.95	8.4	82.4	
L4709-17	MA42Y3	17	1.00	2.00	2.00	100.0	P.E. SAMPLE
L4709-18	MA42Y7	18	1.15	9.75	8.69	87.7	

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$