



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 11/17/2020

OVENTEMP IN Celsius(°C): 107
Time IN: 17:15
In Date: 11/16/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:06
Out Date: 11/17/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB111990

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
L4764-01	WH-PIPE-A	1	1.00	2.00	2.00	100.0	WIPE SAMPLE
L4764-02	WH-PIPE-B	2	1.00	2.00	2.00	100.0	WIPE SAMPLE
L4765-01	CF-OD	3	1.00	2.00	2.00	100.0	OILY DEBRIS
L4765-03	CF-OIL	4	1.00	2.00	2.00	100.0	OIL SAMPLE
L4766-01	PIPE-AUD-31334	5	1.00	2.00	2.00	100.0	WIPE SAMPLE
L4766-02	PIPE-AUD-31335	6	1.00	2.00	2.00	100.0	WIPE SAMPLE
L4766-03	PIPE-AUD-31336	7	1.00	2.00	2.00	100.0	WIPE SAMPLE
L4767-01	NWB-1555	8	1.00	2.00	2.00	100.0	OIL SAMPLE
L4768-01	TR-05-111620	45	1.14	9.59	8.41	86.0	
L4768-02	TR-05-111620-EPH-2	46	1.14	9.94	8.61	84.9	
L4769-01	TRE-1120-1122-COMP	47	1.14	9.73	9.44	96.6	
L4769-07	TRE-1089	48	1.15	9.74	8.75	88.5	
L4769-10	ETGI-359	49	1.16	9.79	9.23	93.5	
L4770-01	FJ-BH-01-0-1	9	1.13	9.9	8.8	87.5	
L4770-02	FJ-BH-01-0-1	10	1.15	9.56	8.96	92.9	
L4770-03	FJ-BH-01-2-3	11	1.16	9.78	8.74	87.9	
L4770-04	FJ-BH-01-2-3	12	1.15	9.51	8.44	87.2	
L4770-05	FJ-BH-01-4-5	13	1.17	9.59	9.33	96.9	
L4770-06	FJ-BH-01-4-5	14	1.16	9.59	9.41	97.9	
L4770-07	FJ-BH-02-0-1	15	1.17	9.6	8.54	87.4	
L4770-08	FJ-BH-02-0-1	16	1.15	9.69	9.06	92.6	
L4770-09	FJ-BH-02-2-3	17	1.16	9.94	8.52	83.8	
L4770-10	FJ-BH-02-2-3	18	1.15	9.72	8.5	85.8	
L4770-11	FJ-BH-02-4-5	19	1.17	9.82	9.53	96.6	
L4770-12	FJ-BH-02-4-5	20	1.15	9.82	9.58	97.2	
L4770-13	FJ-BH-03-0-1	21	1.17	9.82	9.13	92.0	
L4770-14	FJ-BH-03-0-1	22	1.15	9.98	9.17	90.8	
L4770-15	FJ-BH-03-2-3	23	1.16	9.86	9.1	91.3	



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 11/17/2020

OVENTEMP IN Celsius(°C): 107
Time IN: 17:15
In Date: 11/16/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:06
Out Date: 11/17/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB111990

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
L4770-16	FJ-BH-03-2-3	24	1.17	9.88	8.93	89.1	
L4770-17	FJ-BH-03-4-5	25	1.17	9.8	9.44	95.8	
L4770-18	FJ-BH-03-4-5	26	1.15	9.82	9.28	93.8	
L4770-19	FJ-BH-04-0-1	27	1.17	9.95	9.1	90.3	
L4770-20	FJ-BH-04-0-1	28	1.15	9.72	8.79	89.1	
L4770-21	FJ-BH-04-2-3	29	1.15	9.77	8.92	90.1	
L4770-22	FJ-BH-04-2-3	30	1.17	9.75	9.03	91.6	
L4770-23	FJ-BH-04-4-5	31	1.16	9.46	9.14	96.1	
L4770-24	FJ-BH-04-4-5	32	1.15	9.7	9.4	96.5	
L4771-01	CD-BH-01-0-1	33	1.16	9.98	8.52	83.4	
L4771-02	CD-BH-01-0-1	34	1.17	9.77	7.19	70.0	
L4771-03	CD-BH-01-2-3	35	1.17	9.63	8.25	83.7	
L4771-04	CD-BH-01-2-3	36	1.15	9.58	8.17	83.3	
L4771-05	CD-BH-01-4-5	37	1.17	9.6	8.46	86.5	
L4771-06	CD-BH-01-4-5	38	1.16	9.98	8.7	85.5	
L4771-07	CD-BH-02-0-1	39	1.17	9.63	8.12	82.2	
L4771-08	CD-BH-02-0-1	40	1.15	9.78	8.2	81.7	
L4771-09	CD-BH-02-2-3	41	1.16	9.52	8.09	82.9	
L4771-10	CD-BH-02-2-3	42	1.15	9.7	8.22	82.7	
L4771-11	CD-BH-02-4-5	43	1.15	9.86	8.22	81.2	
L4771-12	CD-BH-02-4-5	44	1.15	9.57	8.19	83.6	

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$