



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 11/19/2020

OVENTEMP IN Celsius(°C): 107
Time IN: 17:05
In Date: 11/18/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:10
Out Date: 11/19/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB112025

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
L4778-01	1840	1	1.00	2.00	2.00	100.0	DIESEL
L4778-04	1811	2	1.00	2.00	2.00	100.0	OIL SAMPLE
L4778-06	VNJ-223	3	1.15	9.81	8.97	90.3	
L4778-07	358	4	1.14	9.74	9.04	91.9	
L4785-01	M-6(0-10)	5	1.11	9.65	8.43	85.7	
L4785-02	M-99(0-10)	6	1.13	9.94	8.59	84.7	
L4786-01	KL-SOIL	7	1.15	9.94	8.55	84.2	
L4787-08	BP-VPB178-DUP-20201116	8	1.00	2.00	2.00	100.0	
L4787-09	BP-VPB178-GW-698-700	9	1.00	2.00	2.00	100.0	
L4791-01	EMS-2-0002	10	1.12	9.77	9.22	93.6	
L4791-02	EMS-2-0204	11	1.15	9.61	8.6	88.1	
L4791-03	EMS-2-0406	12	1.12	9.82	9.5	96.3	
L4791-04	EMS-2-0608	13	1.16	9.61	9.21	95.3	
L4791-05	EMS-2-0810	14	1.18	9.58	9.36	97.4	
L4791-06	EMS-3-0002	15	1.15	9.93	9.27	92.5	
L4791-07	EMS-3-0204	16	1.16	9.77	8.73	87.9	
L4791-08	EMS-3-0406	17	1.18	9.56	8.24	84.2	
L4791-09	EMS-3-0608	18	1.1	9.77	8.84	89.3	
L4791-10	EMS-3-0810	19	1.16	9.73	8.87	90.0	
L4791-11	EMS-4-0002	20	1.19	9.86	8.96	89.6	
L4791-12	EMS-4-0204	21	1.1	9.57	8.94	92.6	
L4791-13	EMS-4-0406	22	1.16	9.74	9.33	95.2	
L4791-14	EMS-4-0608	23	1.12	9.67	9.13	93.7	
L4791-15	EMS-4-0810	24	1.15	9.97	9.55	95.2	
L4791-16	304EM-002-0002	25	1.16	9.59	6.61	64.7	
L4791-17	304EM-002-0204	26	1.1	9.86	8.29	82.1	
L4791-18	304EM-004-0002	27	1.13	9.76	7.37	72.3	
L4791-19	FD-01	28	1.16	9.88	8.33	82.2	



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 11/19/2020

OVENTEMP IN Celsius(°C): 107
Time IN: 17:05
In Date: 11/18/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:10
Out Date: 11/19/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB112025

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
L4792-01	304EM-004-0204	29	1.17	9.88	8.24	81.2	
L4792-02	304EM-003-0002	30	1.15	9.95	8.69	85.7	
L4792-03	L4792-02MS	31	1.15	9.95	8.69	85.7	
L4792-04	L4792-02MSD	32	1.15	9.95	8.69	85.7	
L4792-05	304EM-003-0204	33	1.13	9.65	8.36	84.9	
L4792-06	304EM-001-0002	34	1.14	9.62	8.61	88.1	
L4792-07	304EM-001-0204	35	1.15	9.97	8.52	83.6	
L4792-08	250EM-001-0002	36	1.1	9.58	8.43	86.4	
L4792-09	304EM-006-0002	37	1.13	9.83	8.97	90.1	
L4792-10	304EM-006-0204	38	1.1	9.88	8.48	84.1	
L4792-11	250EM-002-0002	39	1.18	9.86	9.04	90.6	
L4792-12	250EM-002-0204	40	1.14	9.96	9.13	90.6	
L4792-13	250EM-003-0002	41	1.14	9.97	9.26	92.0	
L4792-14	250EM-003-0204	42	1.13	9.61	8.96	92.3	
L4792-15	250EM-004-0002	43	1.12	9.81	9.05	91.3	
L4792-16	250EM-004-0204	44	1.16	9.98	9.22	91.4	
L4795-01	SR-L13-L14-SOUTH-1	45	1.13	9.73	8.02	80.1	
L4795-02	SR-L13-L14-SOUTH-2	46	1.16	9.98	7.00	66.2	
L4795-03	SR-L13-L14-NORTH-1	47	1.16	9.97	9.18	91.0	
L4795-04	SR-L13-L14-NORTH-2	48	1.15	9.83	8.64	86.3	
L4795-05	SR-L13-L14-NORTH-3	49	1.14	9.67	8.48	86.0	
L4795-06	SR-L13-L14-NORTH-4	50	1.13	9.57	8.35	85.5	
L4795-07	SR-L13-L14-NORTH-5	51	1.13	9.71	8.49	85.8	
L4795-08	SR-L13-L14-NORTH-6	52	1.13	9.64	8.3	84.3	
L4795-09	SR-L13-L14-NORTH-7	53	1.13	9.59	8.25	84.2	
L4797-01	R-LA-2-5-WC	54	1.13	9.89	7.92	77.5	
L4797-02	R-LA-2-5	55	1.12	9.9	8.09	79.4	
L4802-01	FK-GF-02-064-A	56	1.15	9.51	8.45	87.3	



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 11/19/2020

OVENTEMP IN Celsius(°C): 107
Time IN: 17:05
In Date: 11/18/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:10
Out Date: 11/19/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB112025

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
L4802-02	FK-GF-02-064-A	57	1.18	9.61	8.43	86.0	
L4802-03	FK-GF-02-064-B	58	1.12	9.87	8.85	88.3	
L4802-04	FK-GF-02-064-B	59	1.13	9.62	8.5	86.8	
L4803-01	FK-BPM-01-C5-D5-C5A-D5A	60	1.16	9.61	9.06	93.5	
L4803-02	FK-BPM-01-C5-D5-C5A-D5A	61	1.17	9.77	9.3	94.5	
L4804-01	FK-BPM-04-012-ABCD	62	1.15	9.54	9.00	93.6	
L4804-02	FK-BPM-04-012-ABCD	63	1.13	9.65	9.08	93.3	
L4805-01	FK-BPM-03-017-ABCD	64	1.1	9.76	9.01	91.3	
L4805-02	FK-BPM-03-017-ABCD	65	1.14	9.94	9.25	92.2	
L4806-01	FK-BPM-03-018-ABCD	66	1.13	9.63	8.87	91.1	
L4806-02	FK-BPM-03-018-ABCD	67	1.13	9.59	8.88	91.6	

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$