



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 12/10/2020

OVENTEMP IN Celsius(°C): 107
Time IN: 17:15
In Date: 11/20/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:15
Out Date: 11/21/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB112070

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
L4592-01	EO-01-112020	37	1.14	9.73	9.03	91.9	
L4592-02	EO-01-112020-EPH-2	38	1.16	9.76	8.89	89.9	
L4835-01	OF-3-20201119	1	1.14	9.97	6.31	58.6	
L4835-02	L4835-01MS	2	1.14	9.97	6.31	58.6	
L4835-03	L4835-01MSD	3	1.14	9.97	6.31	58.6	
L4835-04	DUP03-20201119	4	1.17	9.73	5.75	53.5	
L4836-01	SB-01	5	1.14	9.78	8.71	87.6	
L4836-02	SB-02	6	1.12	9.67	8.12	81.9	
L4836-03	SB-03	7	1.16	9.81	8.82	88.6	
L4836-04	SB-04	8	1.14	9.61	8.29	84.4	
L4836-05	SB-05	9	1.13	9.93	9.06	90.1	
L4836-06	SB-06	10	1.16	9.61	7.88	79.5	
L4836-07	SB-07	11	1.15	9.79	8.37	83.6	
L4836-08	SB-08	12	1.18	9.86	9.19	92.3	
L4836-09	SL-02	13	1.15	9.91	8.11	79.5	
L4836-10	SL-01	14	1.17	9.71	7.01	68.4	
L4837-05	SVOC-GPC-BLANK	15	1.00	2.00	2.00	100.0	
L4837-06	PEST-GPC-BLANK	16	1.00	2.00	2.00	100.0	
L4837-07	PEST-GPC-BLANK-SPIKE	17	1.00	2.00	2.00	100.0	
L4837-08	PCB-GPC-BLANK	18	1.00	2.00	2.00	100.0	
L4837-09	PCB-GPC-BLANK-SPIKE	19	1.00	2.00	2.00	100.0	
L4837-10	SVOC-GPC2-BLANK	20	1.00	2.00	2.00	100.0	
L4837-11	PEST-GPC2-BLANK	21	1.00	2.00	2.00	100.0	
L4837-12	PEST-GCP2-BLANK-SPIKE	22	1.00	2.00	2.00	100.0	
L4837-13	PCB-GPC2-BLANK	23	1.00	2.00	2.00	100.0	
L4837-14	PCB-GPC2-BLANK-SPIKE	24	1.00	2.00	2.00	100.0	
L4838-01	KENS-B2-111920-0-2	27	1.17	9.62	8.38	85.3	
L4838-02	KENS-B2-111920-0-2	28	1.16	9.97	8.83	87.1	



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 12/10/2020

OVENTEMP IN Celsius(°C): 107
Time IN: 17:15
In Date: 11/20/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:15
Out Date: 11/21/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB112070

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
L4838-03	KENS-B2-111920-6-8	29	1.18	9.67	8.7	88.6	
L4838-04	KENS-B2-111920-6-8	30	1.14	9.97	8.65	85.1	
L4839-01	B3-9-10	25	1.18	9.98	9.25	91.7	
L4839-02	B5-7-8	26	1.16	9.97	9.39	93.4	
L4840-01	30739	85	1.14	9.9	9.83	99.2	
L4840-02	30740	86	1.18	9.96	9.93	99.7	
L4842-01	11620	87	1.19	9.78	9.33	94.8	
L4842-02	11620-EPH-2	88	1.14	9.75	9.37	95.6	
L4844-01	OR-02-112020	89	1.1	9.66	9.06	93.0	
L4844-02	OR-02-112020-EPH-2	90	1.18	9.62	9.02	92.9	
L4845-01	PH-CON-A	31	1.00	2.00	2.00	100.0	CONCRETE SAMPLE
L4845-02	PH-CON-B	32	1.00	2.00	2.00	100.0	CONCRETE SAMPLE
L4845-03	PH-TP	33	1.17	9.64	8.81	90.2	
L4846-01	RCA-EDISON-3	73	1.00	2.00	2.00	100.0	STONE SAMPLE,100% SOLIDS
L4847-01	EDISON-1A	74	1.00	2.00	2.00	100.0	WOOD SAMPLE
L4847-02	EDISON-2A	75	1.00	2.00	2.00	100.0	wood sample
L4847-03	EDISON-4A	76	1.00	2.00	2.00	100.0	WOOD SAMPLE
L4848-01	TR5- (1-4) -COMP	39	1.19	9.79	9.12	92.2	
L4848-03	TR5-1-BOTTOM	40	1.17	9.69	9.11	93.2	
L4850-01	WOOSPK-038-0002	41	1.14	9.73	8.88	90.1	
L4850-02	WOOSPK-038-0204	42	1.19	9.51	7.5	75.8	
L4850-05	WOOSPK-032-0002	43	1.17	9.54	8.57	88.4	
L4850-06	WOOSPK-032-0204	44	1.15	9.98	6.36	59.0	
L4850-07	FD-04	45	1.14	9.72	6.58	63.4	
L4850-08	WOOSPK-034-0002	46	1.1	9.61	8.87	91.3	
L4850-09	WOOSPK-034-0204	47	1.14	9.85	7.45	72.4	
L4850-12	WOOSPK-036-0002	48	1.15	9.95	9.37	93.4	
L4850-13	L4850-12MS	49	1.15	9.95	9.37	93.4	



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 12/10/2020

OVENTEMP IN Celsius(°C): 107
Time IN: 17:15
In Date: 11/20/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:15
Out Date: 11/21/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB112070

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
L4850-14	L4850-12MSD	50	1.15	9.95	9.37	93.4	
L4850-15	WOOSPK-036-0204	51	1.12	9.72	8.75	88.7	
L4856-01	OUTFALL-SED-RE	34	1.13	9.93	7.74	75.1	
L4856-02	BG-SED-1-RE	35	1.12	9.67	9.06	92.9	
L4856-03	BG-SED-2-RE	36	1.14	9.94	7.15	68.3	
L4857-01	304EM-STOCKPILE	52	1.14	9.94	7.96	77.5	
L4857-02	340EM-STOCKPILE-EAST	53	1.19	9.68	8.82	89.9	
L4857-03	340EM-STOCKPILE-WEST	54	1.1	9.74	8.61	86.9	
L4858-01	PHC-340EM-012	55	1.16	9.66	7.61	75.9	
L4858-02	PHC-340EM-013	56	1.12	9.56	8.1	82.7	
L4863-01	S4-01-A	57	1.14	9.53	8.45	87.1	
L4863-02	S4-02-A	58	1.12	9.66	8.2	82.9	
L4863-03	S4-05-A	59	1.1	9.64	7.92	79.9	
L4863-04	S4-06-A	60	1.19	9.58	7.99	81.0	
L4863-05	S4-01-B	61	1.15	9.5	8.26	85.1	
L4863-06	S4-02-B	62	1.14	9.84	8.14	80.5	
L4863-07	S4-05-B	63	1.15	9.97	8.9	87.9	
L4863-08	S4-06-B	64	1.19	9.79	8.28	82.4	
L4863-09	S4-01-C	65	1.14	9.91	8.7	86.2	
L4863-10	S4-02-C	66	1.12	9.77	7.84	77.7	
L4863-11	S4-03-C	67	1.14	9.56	7.53	75.9	
L4863-12	S4-04-C	68	1.15	9.97	8.32	81.3	
L4863-13	S4-05-C	69	1.13	9.95	8.22	80.4	
L4863-14	S4-06-C	70	1.13	9.67	8.15	82.2	
L4863-15	DUP-5	71	1.12	9.57	8.4	86.2	
L4863-16	DUP-6	72	1.19	9.74	7.94	78.9	
L4868-01	KENS-B3-112020-0-2	79	1.14	9.68	8.56	86.9	
L4868-02	KENS-B3-112020-0-2	80	1.18	9.98	9.12	90.2	



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 12/10/2020

OVENTEMP IN Celsius(°C): 107
Time IN: 17:15
In Date: 11/20/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:15
Out Date: 11/21/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB112070

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
L4868-03	KENS-B3-112020-8-10	81	1.15	9.95	8.85	87.5	
L4868-04	KENS-B3-112020-8-10	82	1.14	9.96	9.11	90.4	
L4868-05	KENS-112020-BD	83	1.13	9.71	8.54	86.4	
L4868-06	KENS-112020-BD	84	1.19	9.98	9.03	89.2	
L4869-01	FCA-1	77	1.00	2.00	2.00	100.0	CAULKING/CERAMIC
L4869-02	CA-1	78	1.00	2.00	2.00	100.0	CAULKING/CERAMIC

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$