



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 11/24/2020

OVENTEMP IN Celsius(°C): 107
Time IN: 17:00
In Date: 11/23/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:05
Out Date: 11/24/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB112097

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
L4588-01	HD-02-112320	23	1.15	9.97	9.14	90.6	
L4588-02	HD-02-112320-EPH-2	24	1.12	9.92	9.2	91.8	
L4597-01	BU-03-112320	32	1.12	9.51	8.57	88.8	
L4597-02	BU-03-112320-EPH-2	33	1.13	9.55	8.44	86.8	
L4597-03	BU-03-112320	34	1.15	9.94	8.99	89.2	
L4871-01	32783	25	1.16	9.72	6.19	58.8	
L4872-01	SU-04-112320	26	1.19	9.62	8.97	92.3	
L4872-02	SU-04-112320-EPH-2	27	1.16	9.88	9.08	90.8	
L4873-01	HR-01-112320	28	1.14	9.84	9.16	92.2	
L4873-02	HR-01-112320-EPH-2	29	1.14	9.96	9.2	91.4	
L4874-01	72-11933	30	1.14	9.77	8.78	88.5	
L4874-02	72-11933-EPH-2	31	1.13	9.94	9.29	92.6	
L4875-01	TR-08-112320	35	1.18	9.63	8.52	86.9	
L4875-02	TR-08-112320-EPH-2	36	1.13	9.95	8.82	87.2	
L4876-01	TP-2A	1	1.15	9.93	7.38	71.0	
L4876-02	TP-2B	2	1.16	9.73	6.45	61.7	
L4876-03	TP-2C	3	1.16	9.76	6.58	63.0	
L4876-04	TP-2D	4	1.13	9.64	6.64	64.7	
L4877-01	P001-CB001-1224-01	5	1.13	9.94	8.61	84.9	
L4877-02	P001-CB002-1230-01	6	1.15	9.98	8.08	78.5	
L4877-03	P001-CB004-2260-01	7	1.15	9.97	9.06	89.7	
L4877-04	P001-CB005-0304-01	8	1.12	9.72	7.32	72.1	
L4877-05	P001-CB005-1824-01	9	1.1	9.76	7.37	72.4	
L4877-06	P001-CB006-1224-01	10	1.15	9.98	8.53	83.6	
L4877-07	P001-CB007-1218-01	11	1.15	9.75	8.25	82.6	
L4877-08	P001-CB008-3036-01	12	1.15	9.95	7.74	74.9	
L4877-09	L4877-08MS	13	1.15	9.95	7.74	74.9	
L4877-10	L4877-08MSD	14	1.15	9.95	7.74	74.9	



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 11/24/2020

OVENTEMP IN Celsius(°C): 107
Time IN: 17:00
In Date: 11/23/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:05
Out Date: 11/24/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB112097

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
L4877-11	P001-CB008-3036-02	15	1.14	9.91	8.27	81.3	
L4877-12	P001-CB012-1224-01	16	1.1	9.56	7.48	75.4	
L4877-13	P001-CB014-1230-01	17	1.15	9.95	7.81	75.7	
L4877-14	P001-CB015-1622-01	18	1.14	9.67	7.49	74.4	
L4877-15	P001-CB015-2228-01	19	1.13	9.64	8.17	82.7	
L4877-16	P001-CB016-1224-01	20	1.14	9.96	8.38	82.1	
L4877-17	P001-CB017-1022-01	21	1.15	9.94	8.66	85.4	
L4877-18	P001-CB018-0924-01	22	1.16	9.98	9.11	90.1	

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$