



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 11/25/2020

OVENTEMP IN Celsius(°C): 107
Time IN: 17:15
In Date: 11/24/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:10
Out Date: 11/25/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB112118

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
L4880-01	KENS-B4-112320-0-2	1	1.12	9.98	8.59	84.3	
L4880-02	KENS-B4-112320-0-2	2	1.13	9.86	8.58	85.3	
L4880-03	L4880-01MS	3	1.12	9.98	8.59	84.3	
L4880-04	L4880-02MS	4	1.13	9.86	8.58	85.3	
L4880-05	L4880-01MSD	5	1.12	9.98	8.59	84.3	
L4880-06	L4880-02MSD	6	1.13	9.86	8.58	85.3	
L4880-07	KENS-B4-112320-6-8	7	1.18	9.69	8.71	88.5	
L4880-08	KENS-B4-112320-6-8	8	1.13	9.59	8.67	89.1	
L4882-01	PAV-B48-112320-0-1	9	1.15	9.55	8.43	86.7	
L4882-03	PAV-B48-112320-0-10	10	1.13	9.81	8.66	86.8	
L4882-05	PAV-B48-112320-10-16.5	11	1.14	9.67	8.65	88.0	
L4882-07	PAV-B48-112320-15.5-16.5	12	1.18	9.91	8.82	87.5	
L4884-01	FS-OD	13	1.00	2.00	2.00	100.0	OILY DEBRIS
L4885-01	OK-03-112420	14	1.15	9.94	9.36	93.4	
L4885-02	OK-03-112420-EPH-2	15	1.16	9.57	8.84	91.3	
L4886-01	BC248138-END-1	16	1.00	2.00	2.00	100.0	CABLE SAMPLE
L4886-02	BC248138-END-2	17	1.00	2.00	2.00	100.0	CABLE SAMPLE
L4886-03	HEC8635-END-1	18	1.00	2.00	2.00	100.0	CABLE SAMPLE
L4886-04	HEC8635-END-2	19	1.00	2.00	2.00	100.0	CABLE SAMPLE
L4889-01	BP-VPB178-GW-768-770	20	1.00	2.00	2.00	100.0	
L4889-02	BP-VPB178-GW-798-800	21	1.00	2.00	2.00	100.0	
L4889-03	BP-VPB178-GW-818-820	22	1.00	2.00	2.00	100.0	
L4890-01	K205-9B-1	23	1.15	9.81	8.52	85.1	
L4890-02	K205-9B-2	24	1.12	9.81	8.34	83.1	
L4891-01	PAV-B3-112420-0-1	25	1.15	9.95	9.14	90.8	
L4891-03	PAV-B3-112420-0-10	26	1.17	9.81	8.62	86.2	
L4891-05	PAV-B3-112420-10-11.5	27	1.19	9.66	8.59	87.4	
L4891-07	PAV-B3-112420-10.5-11.5	28	1.19	9.85	8.8	87.9	



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 11/25/2020

OVENTEMP IN Celsius(°C): 107
Time IN: 17:15
In Date: 11/24/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:10
Out Date: 11/25/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB112118

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
L4892-01	SB-45 (8.0-8.5)	29	1.14	9.94	9.32	93.0	
L4892-04	SB-46 (8.0-8.5)	30	1.14	9.81	8.98	90.4	
L4892-07	SB-47 (8.0-8.5)	31	1.16	9.98	8.48	83.0	
L4892-10	SB-48 (8.0-8.5)	32	1.12	9.77	9.06	91.8	
L4892-13	SB-49 (8.0-8.5)	33	1.16	9.83	8.35	82.9	
L4892-17	SB-50 (8.0-8.5)	34	1.14	9.93	8.72	86.2	
L4892-21	SB-51 (8.0-8.5)	35	1.14	9.97	9.01	89.1	
L4892-25	SB-52 (8.0-8.5)	36	1.17	9.94	8.53	83.9	
L4892-29	DUP-03 (11-23-2020)	37	1.12	9.6	8.05	81.7	
L4893-01	R-LA-36-2-WC	38	1.13	9.67	7.78	77.9	
L4893-02	R-LA-36-2	39	1.14	9.77	7.85	77.8	
L4894-01	R-LA-36-4-WC	40	1.14	9.74	8.11	81.0	
L4894-02	R-LA-36-4	41	1.1	9.91	8.41	83.0	

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$