



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 11/30/2020

OVENTEMP IN Celsius(°C): 107
Time IN: 16:10
In Date: 11/25/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:00
Out Date: 11/26/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB112146

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
L4897-01	88-CHOPIN	1	1.00	2.00	2.00	100.0	OIL SAMPLE
L4897-02	127-CHOPIN	2	1.00	2.00	2.00	100.0	OIL SAMPLE
L4899-01	TP-1	3	1.11	9.64	8.32	84.5	
L4899-02	TP-2	4	1.16	9.85	8.45	83.9	
L4899-03	TP-1-2	5	1.15	9.96	8.85	87.4	
L4900-01	PAV-B2-112520-0-1	6	1.15	9.95	8.93	88.4	
L4900-03	PAV-B2-112520-0-10	7	1.19	9.71	8.65	87.6	
L4900-05	PAV-B2-112520-10-12.75	8	1.16	9.71	8.74	88.7	
L4900-07	PAV-B2-112520-12-12.75	9	1.15	9.94	9.22	91.8	
L4903-01	TP-1	10	1.19	9.91	8.21	80.5	
L4903-02	TP-2	11	1.16	9.83	8.99	90.3	
L4903-03	TP-1-2-COMP	12	1.17	9.66	8.33	84.3	
L4903-05	TP-4	13	1.15	9.78	8.36	83.5	
L4903-06	TP-3	14	1.17	9.98	8.5	83.2	
L4903-07	TP-3-4-COMP	15	1.17	9.9	8.53	84.3	

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$