



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 12/3/2020

OVENTEMP IN Celsius(°C): 107
Time IN: 15:35
In Date: 12/02/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:06
Out Date: 12/03/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB112196

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
L4927-01	MGB1B5	1	1.15	9.6	9.37	97.3	
L4927-02	MGB1B5D	2	1.15	9.6	9.37	97.3	
L4927-03	MGB1B5S	3	1.15	9.6	9.37	97.3	
L4927-04	MGB1B6	4	1.1	9.75	9.55	97.7	
L4927-05	MGB1B7	5	1.17	9.7	9.35	95.9	
L4927-06	MGB1B8	6	1.13	9.57	9.38	97.7	
L4927-07	MGB1B9	7	1.14	9.91	9.71	97.7	
L4927-08	MGB1C0	8	1.13	9.55	9.34	97.5	
L4927-09	MGB1C1	9	1.13	9.92	9.66	97.0	
L4927-10	MGB1C2	10	1.13	9.6	9.48	98.6	
L4927-11	MGB1C3	11	1.14	9.7	9.34	95.8	
L4927-12	MGB1C4	12	1.16	9.71	9.45	97.0	
L4927-13	MGB1C5	13	1.16	9.82	9.77	99.4	
L4927-14	MGB1C6	14	1.00	2.00	2.00	100.0	P.E. SAMPLE

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$