



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 12/7/2020

OVENTEMP IN Celsius(°C): 107
Time IN: 15:30
In Date: 12/04/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M Oven#1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:00
Out Date: 12/05/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-Oven

QC:LB112234

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
L4956-13	X2R72	1	1.00	2.00	2.00	100.0	P.T. SAMPLE
L4956-14	X2R72D	2	1.00	2.00	2.00	100.0	P.T. SAMPLE
L4956-15	X2R72S	3	1.00	2.00	2.00	100.0	P.T. SAMPLE
L4956-16	X2R73	4	1.00	2.00	2.00	100.0	P.T. SAMPLE
L4956-17	X2R73D	5	1.00	2.00	2.00	100.0	P.T. SAMPLE
L4956-18	X2R73S	6	1.00	2.00	2.00	100.0	P.T. SAMPLE
L4956-19	X2R75	7	1.00	2.00	2.00	100.0	P.T. SAMPLE
L4956-20	X2R75D	8	1.00	2.00	2.00	100.0	P.T. SAMPLE
L4956-21	X2R75S	9	1.00	2.00	2.00	100.0	P.T. SAMPLE
L4956-22	X2R76	10	1.00	2.00	2.00	100.0	P.T. SAMPLE
L4956-23	X2R76D	11	1.00	2.00	2.00	100.0	P.T. SAMPLE
L4956-24	X2R76S	12	1.00	2.00	2.00	100.0	P.T. SAMPLE

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$