



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 12/9/2020

OVENTEMP IN Celsius(°C): 107
Time IN: 14:25
In Date: 12/08/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:03
Out Date: 12/09/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB112265

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
L4979-01	GB1F3	1	1.16	9.74	8.85	89.6	
L4979-02	GB1F4	2	1.16	9.82	8.07	79.8	
L4979-03	GB1F5	3	1.15	9.98	9.00	88.9	
L4979-04	GB1F6	4	1.16	9.97	8.78	86.5	
L4979-05	GB1F7	5	1.18	9.85	8.44	83.7	
L4979-06	GB1F8	6	1.13	9.9	8.59	85.1	
L4979-07	GB1F8MS	7	1.13	9.9	8.59	85.1	
L4979-08	GB1F8MSD	8	1.13	9.9	8.59	85.1	
L4979-09	GB1F9	9	1.18	9.66	8.16	82.3	
L4979-10	GB1G0	10	1.11	9.73	8.35	84.0	
L4979-11	GB1G1	11	1.16	9.84	8.51	84.7	
L4979-12	GB1G2	12	1.19	9.97	8.6	84.4	
L4979-13	GB1G3	13	1.00	2.00	2.00	100.0	P.E. SAMPLE
L4979-25	VHBLK02	14	1.00	2.00	2.00	100.0	VHBLK

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$