



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 12/11/2020

OVENTEMP IN Celsius(°C): 107
Time IN: 17:20
In Date: 12/10/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:22
Out Date: 12/11/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB112300

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
L5009-01	KENS-B5-120920-0-2	1	1.15	9.72	8.79	89.1	
L5009-02	KENS-B5-120920-0-8	2	1.15	9.72	8.79	89.1	
L5009-03	KENS-B5-120920-13-15	3	1.15	9.98	9.00	88.9	
L5009-04	KENS-B5-120920-13-15	4	1.15	9.98	9.00	88.9	
L5009-05	KENS-B5-120920-DP	5	1.11	9.94	9.06	90.0	
L5009-06	KENS-B5-120920-DP	6	1.11	9.94	9.06	90.0	
L5010-01	JBD-SB-01-1.0-3.0	7	1.15	9.95	8.56	84.2	
L5010-02	JBD-SB-02-2.0-3.5	8	1.1	9.69	8.28	83.6	
L5010-03	JBD-SB-03-1.0-3.0	9	1.16	9.82	8.22	81.5	
L5010-04	JBD-IB-08-3.0-5.0	10	1.16	9.54	8.19	83.9	
L5010-05	JBD-DUP-120920	11	1.15	9.55	8.2	83.9	
L5013-01	12-10-A	65	1.00	2.00	2.00	100.0	plastic sample
L5013-02	12-10-B	66	1.00	2.00	2.00	100.0	plastic sample
L5013-03	12-10-C	67	1.00	2.00	2.00	100.0	plastic sample
L5013-04	12-10-D	68	1.00	2.00	2.00	100.0	plastic sample
L5013-05	12-10-E	69	1.00	2.00	2.00	100.0	plastic sample
L5013-06	12-10-F	70	1.00	2.00	2.00	100.0	plastic sample
L5013-07	12-10-G	71	1.00	2.00	2.00	100.0	plastic sample
L5013-08	12-10-H	72	1.00	2.00	2.00	100.0	plastic sample
L5013-09	12-10-I	73	1.00	2.00	2.00	100.0	plastic sample
L5013-10	12-10-J	74	1.00	2.00	2.00	100.0	plastic sample
L5013-11	12-10-K	75	1.00	2.00	2.00	100.0	plastic sample
L5013-12	12-10-L	76	1.00	2.00	2.00	100.0	plastic sample
L5013-13	12-10-M	77	1.00	2.00	2.00	100.0	plastic sample
L5014-01	351	49	1.15	9.86	8.95	89.6	
L5014-02	351-EPH-2	50	1.15	9.56	8.76	90.5	
L5015-01	TR-05-121020	51	1.16	9.62	8.46	86.3	
L5015-02	TR-05-121020	52	1.17	9.74	9.1	92.5	



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 12/11/2020

OVENTEMP IN Celsius(°C): 107
Time IN: 17:20
In Date: 12/10/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:22
Out Date: 12/11/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB112300

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
L5015-03	TR-05-121020-EPH-2	53	1.19	9.82	9.26	93.5	
L5016-01	AU-06-121020	54	1.18	9.98	8.94	88.2	
L5016-02	AU-06-121020	55	1.14	9.77	8.77	88.4	
L5016-03	AU-06-121020-EPH-2	56	1.19	9.66	8.65	88.1	
L5017-01	DZA079H-END-1	57	1.00	2.00	2.00	100.0	PILC REEL
L5017-02	DZA079H-END-2	58	1.00	2.00	2.00	100.0	PILC REEL
L5017-03	KMA723H-END-1	59	1.00	2.00	2.00	100.0	PILC REEL
L5017-04	KMA723H-END-2	60	1.00	2.00	2.00	100.0	PILC REEL
L5017-05	GJA307V-END-1	61	1.00	2.00	2.00	100.0	PILC REEL
L5017-06	GJA307V-END-2	62	1.00	2.00	2.00	100.0	PILC REEL
L5018-01	KZZ093L-END-1	45	1.00	2.00	2.00	100.0	PILC REEL
L5018-02	KZZ093L-END-2	46	1.00	2.00	2.00	100.0	PILC REEL
L5018-03	KZ093L-END-1	47	1.00	2.00	2.00	100.0	PILC REEL
L5018-04	KZ093L-END-2	48	1.00	2.00	2.00	100.0	PILC REEL
L5019-01	HIH021D-END-1	63	1.00	2.00	2.00	100.0	PILC REEL
L5019-02	HIH021D-END-2	64	1.00	2.00	2.00	100.0	PILC REEL
L5026-01	SB-53 (4.5-5.0)	24	1.15	9.67	9.06	92.8	
L5026-02	SB-53 (6.5-7.0)	25	1.14	9.75	9.14	92.9	
L5026-03	SB-53 (8.0-8.5)	26	1.16	9.97	9.36	93.1	
L5026-04	SB-54 (4.5-5.0)	27	1.15	9.85	9.23	92.9	
L5026-05	SB-54 (6.5-7.0)	28	1.15	9.93	9.5	95.1	
L5026-06	SB-54 (8.0-8.5)	29	1.15	9.89	9.18	91.9	
L5026-07	SB-55 (4.5-5.0)	30	1.12	9.83	9.14	92.1	
L5026-08	SB-55 (6.5-7.0)	31	1.14	9.98	9.1	90.0	
L5026-09	SB-55 (8.0-8.5)	32	1.13	9.72	8.12	81.4	
L5026-16	SB-58 (4.5-5.0)	33	1.13	9.73	8.76	88.7	
L5026-17	SB-58 (6.5-7.0)	34	1.15	9.95	9.26	92.2	
L5026-18	SB-58 (8.0-8.5)	35	1.14	9.9	9.36	93.8	



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 12/11/2020

OVENTEMP IN Celsius(°C): 107
Time IN: 17:20
In Date: 12/10/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:22
Out Date: 12/11/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB112300

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
L5026-29	SB-62 (4.5-5.0)	36	1.13	9.81	8.88	89.3	
L5026-30	SB-62 (6.5-7.0)	37	1.14	9.83	8.87	89.0	
L5026-31	SB-62 (8.0-8.5)	38	1.14	9.95	9.51	95.0	
L5026-32	SB-63 (4.5-5.0)	39	1.14	9.98	9.44	93.9	
L5026-33	SB-63 (6.5-7.0)	40	1.1	9.88	9.12	91.3	
L5026-34	SB-63 (8.0-8.5)	41	1.14	9.8	9.02	91.0	
L5026-35	SB-64 (4.5-5.0)	42	1.16	9.68	8.89	90.7	
L5026-36	SB-64 (6.5-7.0)	43	1.17	9.69	8.83	89.9	
L5026-37	SB-64 (8.0-8.5)	44	1.14	9.96	9.19	91.3	
L5033-01	CS-43	12	1.13	9.91	9.6	96.5	
L5033-02	CS-44	13	1.16	9.66	9.36	96.5	
L5033-03	CS-45	14	1.11	9.68	9.34	96.0	
L5033-04	CS-46	15	1.13	9.76	9.44	96.3	
L5033-05	CS-47	16	1.19	9.9	9.6	96.6	
L5033-06	CS-48	17	1.15	9.66	9.32	96.0	
L5033-06D	CS-48DUP	18	1.17	9.69	9.35	96.0	
L5035-01	PE-7	19	1.13	9.54	7.95	81.1	
L5035-02	PE-1A	20	1.18	9.59	8.41	86.0	
L5035-03	PE-6	21	1.15	9.98	8.8	86.6	
L5035-04	PE-8	22	1.15	9.95	9.22	91.7	
L5035-05	DUP	23	1.13	9.62	7.88	79.5	

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$