



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 12/16/2020

OVENTEMP IN Celsius(°C): 107
Time IN: 15:25
In Date: 12/15/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:03
Out Date: 12/16/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB112391

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
L5071-01	A54J1	1	1.15	9.65	7.95	80.0	
L5071-02	A54J2	2	1.18	9.83	7.64	74.7	
L5071-03	A54J3	3	1.15	9.98	8.15	79.3	
L5071-04	A54J4	4	1.13	9.51	7.69	78.3	
L5071-05	A54J5	5	1.12	9.64	7.95	80.2	
L5071-06	A54J6	6	1.13	9.93	7.95	77.5	
L5071-07	A54J7	7	1.14	9.97	8.19	79.8	
L5071-08	A54J8	8	1.13	9.74	7.99	79.7	
L5071-09	A54J9	9	1.19	9.58	8.12	82.6	
L5071-10	A54K0	10	1.15	9.95	8.66	85.3	
L5071-11	A54K1	11	1.1	9.86	8.87	88.7	
L5071-13	A54K5	12	1.00	2.00	2.00	100.0	P.E. SAMPLE
L5071-14	A54K7	13	1.1	9.68	8.97	91.7	
L5071-15	A54K8	14	1.12	9.7	8.74	88.8	
L5071-16	A54K9	15	1.1	9.88	9.16	91.8	
L5071-17	A54L0	16	1.13	9.58	8.63	88.8	
L5071-18	A54L1	17	1.16	9.66	8.57	87.2	
L5071-19	A54L1MS	18	1.16	9.66	8.57	87.2	
L5071-20	A54L1MSD	19	1.16	9.66	8.57	87.2	

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$