



PERCENT SOLID

Supervisor: apatel
Analyst: JIGNESH
Date: 12/17/2020

OVENTEMP IN Celsius(°C): 107
Time IN: 16:30
In Date: 12/16/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:15
Out Date: 12/17/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB112401

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
L5120-01	K650-1A	1	1.15	9.64	7.08	69.8	
L5120-02	K650-1B	2	1.16	9.54	7.63	77.2	
L5120-03	K650-1C	3	1.18	9.66	7.58	75.5	
L5120-04	K650-2B	4	1.13	9.66	7.94	79.8	
L5120-05	K650-2C	5	1.16	9.81	7.4	72.1	
L5120-06	K650-3B	6	1.16	9.97	7.61	73.2	
L5120-07	K650-3C	7	1.15	9.92	7.65	74.1	
L5120-08	K650-4B	8	1.13	9.96	7.61	73.4	
L5120-09	K650-4C	9	1.15	9.96	7.65	73.8	
L5120-10	K650-5A	10	1.14	9.98	7.91	76.6	
L5120-11	K650-5B	11	1.13	9.94	7.82	75.9	
L5120-12	K650-5C	12	1.13	9.91	7.9	77.1	
L5120-13	K650-6A	13	1.14	9.5	7.07	70.9	
L5120-14	K650-6B	14	1.13	9.85	7.37	71.6	
L5120-15	K650-6C	15	1.16	9.83	7.79	76.5	
L5120-16	K650-7A	16	1.13	9.75	7.5	73.9	
L5120-17	K650-7B	17	1.16	9.62	7.25	72.0	
L5120-18	K650-7C	18	1.14	9.82	7.68	75.3	
L5120-19	K650-8C	19	1.17	9.66	7.87	78.9	
L5120-20	K650-9C	20	1.16	9.96	8.09	78.7	
L5121-01	K650-10A	21	1.14	9.89	7.8	76.1	
L5121-01D	K650-10ADUP	22	1.15	9.91	7.82	76.1	
L5121-04	K650-11A	23	1.16	9.67	7.75	77.4	
L5121-07	K650-12A	24	1.14	9.67	7.58	75.5	
L5121-10	K650-13A	25	1.13	9.5	7.25	73.1	
L5121-13	K650-14A	26	1.16	9.81	7.28	70.8	
L5121-16	K650-15A	27	1.16	9.9	7.96	77.8	
L5122-01	K650-16A	28	1.14	9.89	8.00	78.4	



PERCENT SOLID

Supervisor: apatel
Analyst: JIGNESH
Date: 12/17/2020

OVENTEMP IN Celsius(°C): 107
Time IN: 16:30
In Date: 12/16/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:15
Out Date: 12/17/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB112401

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
L5122-04	K650-17A	29	1.17	9.85	7.43	72.1	
L5122-07	K650-18A	30	1.15	9.63	7.71	77.4	
L5122-10	K650-19A	31	1.1	9.86	7.76	76.0	
L5122-13	K650-20A	32	1.13	9.64	8.66	88.5	
L5122-15	K650-21A	33	1.14	9.81	8.24	81.9	
L5122-17	K650-DUP1	34	1.16	9.5	7.4	74.8	
L5122-18	K650-DUP2	35	1.18	9.64	7.57	75.5	
L5122-19	K650-DUP3	36	1.15	9.96	8.07	78.5	
L5122-22	K650-DUP6	37	1.12	9.55	8.62	89.0	
L5122-22D	K650-DUP6DUP	38	1.13	9.57	8.64	89.0	
L5125-01	DETENTION-LINE-1	39	1.16	9.71	8.83	89.7	
L5125-02	DETENTION-LINE-1	40	1.13	9.81	8.6	86.1	
L5125-03	DETENTION-LINE-2	41	1.15	9.95	9.24	91.9	
L5125-04	DETENTION-LINE-2	42	1.17	9.9	8.86	88.1	
L5125-05	DETENTION-LINE-3	43	1.16	9.56	8.18	83.6	
L5125-06	DETENTION-LINE-3	44	1.14	9.98	9.03	89.3	

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$