



PERCENT SOLID

Supervisor: apatel
Analyst: JIGNESH
Date: 12/21/2020

OVENTEMP IN Celsius(°C): 107
Time IN: 15:40
In Date: 12/18/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:05
Out Date: 12/19/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB112440

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
L5149-01	CB897	1	1.14	9.7	9.58	98.6	
L5149-02	CB898	2	1.15	9.56	9.43	98.5	
L5149-03	CB898MS	3	1.15	9.56	9.43	98.5	
L5149-04	CB898MSD	4	1.15	9.56	9.43	98.5	
L5149-05	CB899	5	1.15	9.56	9.49	99.2	
L5149-06	CB8A0	6	1.16	9.89	9.72	98.1	
L5149-07	CB8A1	7	1.12	9.52	9.42	98.8	
L5149-08	CB8A2	8	1.18	9.65	9.63	99.8	

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$

12/18/2020

WORKLIST(Hardcopy Internal Chain)

WorkList Name : %1-L5149

WorkList ID : 145082

Department : Wet-Chemistry

Date : 12-18-2020 15:13:46

Due Date	Matrix	Sample	Test	Preservative	Customer	Storage Location	Customer Sample	Collect Date	Method
11/20/2020	Solid	L5149-01	Percent Solids	Cool 4 deg C	USEP04	B41	CB897	11/10/2020	Chemtech -SO
12/17/2020	Solid	L5149-02	Percent Solids	Cool 4 deg C	USEP04	B41	CB898	12/07/2020	Chemtech -SO
12/17/2020	Solid	L5149-03	Percent Solids	Cool 4 deg C	USEP04	B41	CB898MS	12/07/2020	Chemtech -SO
12/17/2020	Solid	L5149-04	Percent Solids	Cool 4 deg C	USEP04	B41	CB898MSD	12/07/2020	Chemtech -SO
12/17/2020	Solid	L5149-05	Percent Solids	Cool 4 deg C	USEP04	B41	CB899	12/07/2020	Chemtech -SO
12/17/2020	Solid	L5149-06	Percent Solids	Cool 4 deg C	USEP04	B41	CB8A0	12/07/2020	Chemtech -SO
12/17/2020	Solid	L5149-07	Percent Solids	Cool 4 deg C	USEP04	B41	CB8A1	12/07/2020	Chemtech -SO
12/25/2020	Solid	L5149-08	Percent Solids	Cool 4 deg C	USEP04	B41	CB8A2	12/15/2020	Chemtech -SO

Date/Time 12-18-2020 15:20
 Received by: J. C. Wet-Chem
 Relinquished by: J. C. Wet-Chem

Date/Time 12-18-2020 16:00
 Received by: J. C. Wet-Chem
 Relinquished by: J. C. Wet-Chem