



PERCENT SOLID

Supervisor: apatel
Analyst: JIGNESH
Date: 12/23/2020

OVENTEMP IN Celsius(°C): 107
Time IN: 17:50
In Date: 12/22/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:25
Out Date: 12/23/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB112487

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
L5166-01	S6-01-A	1	1.15	9.73	7.41	73.0	
L5166-02	S6-02-A	2	1.18	9.78	8.03	79.7	
L5166-03	S6-03-A	3	1.16	9.7	7.63	75.8	
L5166-04	S6-04-A	4	1.15	9.94	8.37	82.1	
L5166-05	S6-05-A	5	1.16	9.63	8.24	83.6	
L5166-06	S6-06-A	6	1.16	9.77	7.73	76.3	
L5166-07	S6-07-A	7	1.14	9.95	8.37	82.1	
L5166-08	S6-08-A	8	1.13	9.5	7.88	80.6	
L5166-09	S6-09-A	9	1.15	9.96	5.65	51.1	
L5166-10	S6-10-A	10	1.14	9.91	8.19	80.4	
L5166-11	S6-11-A	11	1.12	9.74	7.92	78.9	
L5166-12	S6-06-B	12	1.15	9.95	6.99	66.4	
L5166-13	S6-07-B	13	1.14	9.95	8.2	80.1	
L5166-14	S6-08-B	14	1.19	9.82	7.96	78.4	
L5166-15	S6-09-B	15	1.16	9.8	7.77	76.5	
L5166-16	S6-10-B	16	1.14	9.82	7.61	74.5	
L5166-17	S6-11-B	17	1.16	9.72	7.71	76.5	
L5166-18	S6-06-C	18	1.18	9.56	7.56	76.1	
L5166-19	S6-07-C	19	1.14	9.84	8.4	83.4	
L5166-20	S6-08-C	20	1.16	9.6	7.89	79.7	
L5166-20D	S6-08-CDUP	21	1.19	9.62	7.91	79.7	
L5167-01	S7-01-A	22	1.18	9.68	6.81	66.2	
L5167-02	S7-02-A	23	1.16	9.87	7.00	67.0	
L5167-03	S7-03-A	24	1.16	9.74	7.39	72.6	
L5167-04	S8-01-A	25	1.17	9.7	7.4	73.0	
L5167-05	S8-02-A	26	1.14	9.91	8.4	82.8	
L5167-06	S8-03-A	27	1.16	9.76	8.15	81.3	
L5167-07	S8-04-A	28	1.17	9.71	8.21	82.4	



PERCENT SOLID

Supervisor: apatel
Analyst: JIGNESH
Date: 12/23/2020

OVENTEMP IN Celsius(°C): 107
Time IN: 17:50
In Date: 12/22/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:25
Out Date: 12/23/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB112487

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
L5167-08	S8-01-B	29	1.15	9.98	8.45	82.7	
L5167-09	S8-02-B	30	1.14	9.5	7.84	80.1	
L5167-10	S8-03-B	31	1.12	9.6	8.11	82.4	
L5167-11	S8-04-B	32	1.16	9.77	8.3	82.9	
L5167-12	S8-01-C	33	1.15	9.97	7.61	73.2	
L5167-13	S8-02-C	34	1.16	9.68	8.18	82.4	
L5167-14	S8-03-C	35	1.14	9.98	8.35	81.6	
L5167-15	S8-04-C	36	1.16	9.75	8.12	81.0	
L5167-16	S8-05-C	37	1.14	9.57	8.03	81.7	
L5167-17	S8-06-C	38	1.12	9.66	7.03	69.2	
L5167-18	DUP-13	39	1.15	9.97	7.04	66.8	
L5167-19	DUP-14	40	1.14	9.74	7.91	78.7	
L5167-19D	DUP-14DUP	41	1.16	9.76	7.93	78.7	
L5167-20	S9-01-C	42	1.12	9.72	7.98	79.8	
L5167-21	DUP-15	43	1.13	9.61	8.14	82.7	
L5168-01	S13-01-C	44	1.14	9.95	6.78	64.0	
L5168-02	S13-02-C	45	1.15	9.98	7.5	71.9	
L5168-03	S13-03-C	46	1.14	9.94	7.69	74.4	
L5168-04	S13-04-C	47	1.15	9.5	4.05	34.7	
L5168-05	S14-01-A	48	1.17	9.82	8.04	79.4	
L5168-06	S14-02-A	49	1.19	9.81	8.00	79.0	
L5168-07	DUP-27	50	1.1	9.94	8.33	81.8	
L5168-08	S15-01-A	51	1.18	9.58	7.72	77.9	
L5168-09	S15-02-A	52	1.15	9.96	8.12	79.1	
L5168-10	S16-01-A	53	1.18	9.96	5.8	52.6	
L5168-11	S16-02-A	54	1.1	9.92	6.52	61.5	
L5168-12	DUP-30	55	1.14	9.6	7.42	74.2	
L5168-13	S17-01-A	56	1.18	9.58	7.52	75.5	



PERCENT SOLID

Supervisor: apatel
Analyst: JIGNESH
Date: 12/23/2020

OVENTEMP IN Celsius(°C): 107
Time IN: 17:50
In Date: 12/22/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:25
Out Date: 12/23/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB112487

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
L5168-14	S17-02-A	57	1.12	9.8	7.76	76.5	
L5168-15	S17-03-A	58	1.18	9.69	8.51	86.1	
L5168-16	S17-04-A	59	1.16	9.62	8.04	81.3	
L5168-17	S17-05-A	60	1.15	9.69	8.33	84.1	
L5168-17D	S17-05-ADUP	61	1.17	9.71	8.35	84.1	
L5168-18	S17-06-A	62	1.15	9.97	9.04	89.5	
L5168-19	S17-07-A	63	1.13	9.62	9.00	92.7	
L5169-01	S18-01-A	64	1.14	9.68	7.35	72.7	
L5169-02	S18-02-A	65	1.18	9.83	8.11	80.1	
L5169-03	S18-03-A	66	1.14	9.84	7.75	76.0	
L5169-04	S18-04-A	67	1.19	9.62	7.72	77.5	
L5169-05	S18-05-A	68	1.11	9.72	7.27	71.5	
L5169-06	S18-06-A	69	1.15	9.67	7.82	78.3	
L5169-07	S18-01-B	70	1.12	9.55	7.39	74.4	
L5169-08	S18-02-B	71	1.19	9.63	7.77	78.0	
L5169-09	S18-03-B	72	1.18	9.77	6.71	64.4	
L5169-10	S18-04-B	73	1.19	9.75	8.99	91.1	
L5169-11	S18-05-B	74	1.11	9.88	7.74	75.6	
L5169-12	S18-06-B	75	1.16	9.54	6.86	68.0	
L5169-13	S18-01-C	76	1.19	9.67	7.59	75.5	
L5169-14	S18-02-C	77	1.15	9.94	7.21	68.9	
L5169-15	S18-03-C	78	1.13	9.55	7.49	75.5	
L5169-16	S18-04-C	79	1.18	9.64	8.71	89.0	
L5169-17	S18-05-C	80	1.16	9.68	7.54	74.9	
L5169-17D	S18-05-CDUP	81	1.19	9.72	7.58	74.9	
L5169-18	S18-06-C	82	1.14	9.56	6.02	58.0	
L5169-19	S18-07-C	83	1.19	9.78	5.12	45.8	
L5169-20	S18-08-C	84	1.17	9.7	7.09	69.4	



PERCENT SOLID

Supervisor: apatel
Analyst: JIGNESH
Date: 12/23/2020

OVENTEMP IN Celsius(°C): 107
Time IN: 17:50
In Date: 12/22/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:25
Out Date: 12/23/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB112487

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
L5169-21	S18-09-C	85	1.16	9.97	7.4	70.8	
L5169-22	DUP-33	86	1.14	9.95	7.47	71.9	
L5170-01	S19-01-A	94	1.12	9.87	8.06	79.3	
L5170-02	S19-02-A	95	1.14	9.93	8.13	79.5	
L5170-03	S19-03-A	96	1.14	9.71	7.96	79.6	
L5170-04	S19-04-A	97	1.15	9.93	8.15	79.7	
L5170-05	S19-05-A	98	1.18	9.5	7.67	78.0	
L5170-06	S19-06-A	99	1.18	9.97	8.6	84.4	
L5170-07	S19-07-A	100	1.17	9.65	8.37	84.9	
L5170-07D	S19-07-ADUP	101	1.19	9.69	8.41	84.9	
L5170-08	DUP-36	102	1.19	9.7	8.39	84.6	
L5170-09	DUP-18	103	1.13	9.69	9.3	95.4	
L5170-10	DUP-19	104	1.14	9.97	8.93	88.2	
L5170-11	S11-01-A	105	1.18	9.81	8.73	87.5	
L5170-12	S11-02-A	106	1.18	9.93	8.85	87.7	
L5170-13	S11-03-A	107	1.19	9.88	8.94	89.2	
L5170-14	S11-04-A	108	1.12	9.98	8.99	88.8	
L5170-15	S11-05-A	109	1.17	9.71	8.78	89.1	
L5170-16	S11-06-A	110	1.18	9.94	8.25	80.7	
L5170-17	S11-07-A	111	1.16	9.84	8.36	82.9	
L5170-18	S11-08-A	112	1.13	9.6	7.91	80.0	
L5170-19	S11-09-A	113	1.18	9.97	8.68	85.3	
L5170-20	DUP-20	114	1.13	9.66	8.15	82.3	
L5186-01	SP-138KV-YARD	115	1.14	9.98	9.27	92.0	
L5186-02	SP-138KV-YARD-EPH-2	116	1.13	9.94	8.72	86.2	
L5188-01	NB-08-122220	117	1.17	9.71	8.45	85.2	
L5188-02	NB-08-122220-EPH-2	118	1.12	9.53	8.33	85.7	
L5189-01	22	123	1.14	9.53	7.43	75.0	



PERCENT SOLID

Supervisor: apatel
Analyst: JIGNESH
Date: 12/23/2020

OVENTEMP IN Celsius(°C): 107
Time IN: 17:50
In Date: 12/22/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:25
Out Date: 12/23/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB112487

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
L5189-05	22-EPH-2	124	1.19	9.79	7.78	76.6	
L5189-06	09	125	1.17	9.6	9.11	94.2	
L5189-10	09-EPH-2	126	1.19	9.75	9.06	91.9	
L5189-11	18	127	1.19	9.94	8.92	88.3	
L5189-15	18-EPH-2	128	1.19	9.84	8.63	86.0	
L5189-16	08	129	1.19	9.82	7.65	74.9	
L5189-20	08-EPH-2	130	1.17	9.63	7.48	74.6	
L5190-01	OK-01-122220	119	1.15	9.96	9.26	92.1	
L5190-02	OK-01-122220-EPH-2	120	1.17	9.64	8.88	91.0	
L5192-01	CL-02-122220	121	1.14	9.92	8.83	87.6	
L5192-02	CL-02-122220-EPH-2	122	1.12	9.56	7.91	80.5	
L5193-01	LLEXSITUTV-021-001-SO	87	1.13	9.51	9.5	99.9	
L5193-02	L5193-01MS	88	1.13	9.51	9.5	99.9	
L5193-03	L5193-01MSD	89	1.13	9.51	9.5	99.9	
L5193-04	LLEXSITUTV-022-001-SO	90	1.19	9.66	9.63	99.6	
L5193-05	LLEXSITUTV-022-002-SO	91	1.19	9.74	9.71	99.6	
L5193-06	LLEXSITUTV-023-001-SO	92	1.19	9.61	9.56	99.4	
L5193-07	LLEXSITUTV-024-001-SO	93	1.18	9.67	9.6	99.2	

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$