



PERCENT SOLID

Supervisor: apatel
Analyst: JIGNESH
Date: 12/31/2020

OVENTEMP IN Celsius(°C): 107
Time IN: 16:00
In Date: 12/30/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 07:45
Out Date: 12/31/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB112573

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
L5214-03	MDL-SOIL-QT1-2020-01	1	1.00	2.00	2.00	100.0	
L5214-04	MDL-SOIL-QT1-2020-02	2	1.00	2.00	2.00	100.0	
L5214-05	MDL-MED-SOIL-QT1-2020-01	3	1.00	2.00	2.00	100.0	
L5214-06	MDL-MED-SOIL-QT1-2020-02	4	1.00	2.00	2.00	100.0	
L5242-01	SAMPLE-1 (494+26)	5	1.18	9.78	9.37	95.2	
L5242-02	SAMPLE-1 (494+26) VOC-GRAB	6	1.13	5.59	5.41	96.0	
L5242-03	SAMPLE-2 (484+26)	7	1.17	9.85	9.13	91.7	
L5242-04	SAMPLE-2 (484+26) VOC-GRAB	8	1.12	5.54	5.2	92.3	
L5242-05	SAMPLE-3 (474+26)	9	1.15	9.95	9.41	93.9	
L5242-06	SAMPLE-3 (474+26) VOC-GRAB	10	1.12	5.78	5.23	88.2	
L5242-07	SAMPLE-4 (464+26)	11	1.13	9.9	9.33	93.5	
L5242-08	SAMPLE-4 (464+26) VOC-GRAB	12	1.14	6.12	5.77	93.0	
L5242-09	SAMPLE-5 (454+26)	13	1.16	9.77	9.21	93.5	
L5242-10	SAMPLE-5 (454+26) VOC-GRAB	14	1.12	6.55	6.19	93.4	
L5242-11	SAMPLE-6 (444+26)	15	1.19	9.78	9.34	94.9	
L5242-11D	SAMPLE-6 (444+26) DUP	16	1.17	9.8	9.36	94.9	
L5242-12	SAMPLE-6 (444+26) VOC-GRAB	17	1.13	5.43	5.23	95.3	
L5242-13	SAMPLE-7 (434+26)	18	1.11	9.74	9.19	93.6	
L5242-14	SAMPLE-7 (434+26) VOC-GRAB	19	1.18	5.71	5.42	93.6	
L5244-01	HR-01-123020	20	1.12	9.83	9.2	92.8	
L5244-02	HR-01-123020-EPH-2	21	1.13	9.72	8.86	90.0	
L5246-01	1	22	1.12	9.54	8.88	92.2	
L5246-02	2	23	1.19	9.57	8.82	91.1	



PERCENT SOLID

Supervisor: apatel
Analyst: JIGNESH
Date: 12/31/2020

OVENTEMP IN Celsius(°C): 107
Time IN: 16:00
In Date: 12/30/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 07:45
Out Date: 12/31/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB112573

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
I5246-03	3	24	1.15	9.96	9.44	94.1	
I5246-04	4	25	1.1	9.61	8.97	92.5	
I5246-05	5	26	1.16	9.93	9.24	92.1	
I5246-06	6	27	1.19	9.83	9.15	92.1	

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$