



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 1/6/2021

OVENTEMP IN Celsius(°C): 107
Time IN: 16:15
In Date: 01/05/2021
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:00
Out Date: 01/06/2021
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB112632

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
M1008-01	31375	23	1.00	2.00	2.00	100.0	OILY DEBRIS
M1008-02	31376	24	1.00	2.00	2.00	100.0	OIL SAMPLE
M1008-03	31377	25	1.00	2.00	2.00	100.0	OIL SAMPLE
M1008-04	31378	26	1.00	2.00	2.00	100.0	OIL SAMPLE
M1008-05	31379	27	1.00	2.00	2.00	100.0	OIL SAMPLE
M1008-06	31380	28	1.00	2.00	2.00	100.0	OIL SAMPLE
M1008-07	31381	29	1.00	2.00	2.00	100.0	OIL SAMPLE
M1008-08	31382	30	1.00	2.00	2.00	100.0	OIL SAMPLE
M1008-09	31383	31	1.00	2.00	2.00	100.0	OIL SAMPLE
M1008-10	31384	32	1.00	2.00	2.00	100.0	OIL SAMPLE
M1008-11	31788	33	1.00	2.00	2.00	100.0	OIL SAMPLE
M1008-12	31789	34	1.00	2.00	2.00	100.0	OIL SAMPLE
M1008-13	31790	35	1.00	2.00	2.00	100.0	OIL SAMPLE
M1008-14	31791	36	1.00	2.00	2.00	100.0	OIL SAMPLE
M1008-15	31792	37	1.00	2.00	2.00	100.0	OIL SAMPLE
M1008-16	31793	38	1.00	2.00	2.00	100.0	OILY DEBRIS
M1008-17	375-393	39	1.00	2.00	2.00	100.0	OILY DEBRIS
M1008-19	376-384	40	1.00	2.00	2.00	100.0	OIL SAMPLE
M1008-21	388-392	41	1.00	2.00	2.00	100.0	OIL SAMPLE
M1010-01	MT	1	1.00	2.00	2.00	100.0	OIL SAMPLE
M1012-01	KKA091Y-END-1	2	1.00	2.00	2.00	100.0	PILC REEL
M1012-02	KKA091Y-END-2	3	1.00	2.00	2.00	100.0	PILC REEL
M1012-03	DBA351R-END-1	4	1.00	2.00	2.00	100.0	PILC REEL
M1012-04	DBA351R-END-2	5	1.00	2.00	2.00	100.0	PILC REEL
M1013-01	FDA7010-END-1	6	1.00	2.00	2.00	100.0	PILC REEL
M1013-02	FDA7010-END-2	7	1.00	2.00	2.00	100.0	PILC REEL
M1013-03	HEF198S-END-1	8	1.00	2.00	2.00	100.0	PILC REEL
M1013-04	HEF198S-END-2	9	1.00	2.00	2.00	100.0	PILC REEL



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 1/6/2021

OVENTEMP IN Celsius(°C): 107
Time IN: 16:15
In Date: 01/05/2021
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:00
Out Date: 01/06/2021
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB112632

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
M1014-01	HEE293P-END-1	10	1.00	2.00	2.00	100.0	PILC REEL
M1014-02	HEE293P-END-2	11	1.00	2.00	2.00	100.0	PILC REEL
M1014-03	BC226700LOZ-FULL-END-1	12	1.00	2.00	2.00	100.0	PILC REEL
M1014-04	BC226716LOZ-LITTLE-PIECE	13	1.00	2.00	2.00	100.0	PILC REEL
M1014-05	BC226716LOZ-END-1	14	1.00	2.00	2.00	100.0	PILC REEL
M1014-06	BC226716LOZ-END-2	15	1.00	2.00	2.00	100.0	PILC REEL
M1014-07	94021-END-1-1	16	1.00	2.00	2.00	100.0	PILC REEL
M1014-08	94021-END-1-2	17	1.00	2.00	2.00	100.0	PILC REEL
M1014-09	94021-END-2-1	18	1.00	2.00	2.00	100.0	PILC REEL
M1014-10	94021-END-2-2	19	1.00	2.00	2.00	100.0	PILC REEL
M1015-01	NB-01-010521	20	1.13	9.86	9.15	91.9	
M1015-02	NB-01-010521-EPH-2	21	1.15	9.78	8.79	88.5	
M1015-02D	NB-01-010521-EPH-2DUP	22	1.16	9.8	8.81	88.5	

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$