



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 1/8/2021

OVENTEMP IN Celsius(°C): 107
Time IN: 15:03
In Date: 01/07/2021
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:03
Out Date: 01/08/2021
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB112673

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
M1051-13	X2W77	1	1.00	2.00	2.00	100.0	P.E.SAMPLE
M1051-14	X2W77D	2	1.00	2.00	2.00	100.0	P.E.SAMPLE
M1051-15	X2W77S	3	1.00	2.00	2.00	100.0	P.E.SAMPLE
M1051-16	X2W78	4	1.00	2.00	2.00	100.0	P.E.SAMPLE
M1051-17	X2W78D	5	1.00	2.00	2.00	100.0	P.E.SAMPLE
M1051-18	X2W78S	6	1.00	2.00	2.00	100.0	P.E.SAMPLE
M1051-19	X2W79	7	1.00	2.00	2.00	100.0	P.E.SAMPLE
M1051-20	X2W79D	8	1.00	2.00	2.00	100.0	P.E.SAMPLE
M1051-21	X2W79S	9	1.00	2.00	2.00	100.0	P.E.SAMPLE
M1051-22	X2W80	10	1.00	2.00	2.00	100.0	P.E.SAMPLE
M1051-23	X2W80D	11	1.00	2.00	2.00	100.0	P.E.SAMPLE
M1051-24	X2W80S	12	1.00	2.00	2.00	100.0	P.E.SAMPLE

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$