

**PERCENT SOLID**

Supervisor: apatel
Analyst: jignesh
Date: 1/15/2021

OVENTEMP IN Celsius(°C): 107
Time IN: 17:15
In Date: 01/14/2021
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:13
Out Date: 01/15/2021
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB112785

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
M1100-02	TOPSOIL-RBF-ASHCROFT	1	1.17	9.68	7.72	77.0	
M1100-03	TOPSOIL-RBF-DARK	2	1.16	9.94	7.98	77.7	
M1100-04	TOPSOIL-RBF-LIGHT	3	1.12	9.84	9.06	91.1	
M1139-01	PL-01-011421	4	1.15	9.54	8.72	90.2	
M1139-01D	PL-01-011421DUP	5	1.16	9.56	8.74	90.2	
M1139-02	PL-01-011421-EPH-2	6	1.16	9.64	8.62	88.0	
M1140-01	SU-01-011421	7	1.18	9.97	9.06	89.6	
M1140-02	SU-01-011421-EPH-2	8	1.12	9.71	8.97	91.4	
M1141-01	RBR-200034	10	1.15	9.59	8.55	87.7	
M1141-02	RBR-200047	11	1.14	9.97	8.78	86.5	
M1141-03	72-12016	12	1.15	9.93	9.00	89.4	
M1141-04	VNJ-202	13	1.2	9.74	8.18	81.7	
M1141-05	ETGI-360	14	1.16	9.97	9.03	89.3	
M1141-06	RBR-200036	15	1.16	9.75	8.66	87.3	
M1141-07	RBR-200044	16	1.14	9.9	8.87	88.2	
M1141-08	R-4225	17	1.1	9.94	9.06	90.0	
M1141-09	RBR-200029	18	1.15	9.96	9.4	93.6	
M1141-10	1052	19	1.17	9.96	9.52	95.0	
M1141-11	RBR-200054	20	1.17	9.93	8.97	89.0	
M1141-13	RBR-200030	21	1.13	9.93	9.12	90.8	
M1141-14	VNJ-240	22	1.19	9.98	9.09	89.9	
M1141-15	ETGI-340	23	1.13	9.92	9.18	91.6	
M1141-16	RBR-200031	24	1.13	9.93	8.76	86.7	
M1141-17	RBR-200018	25	1.16	9.96	8.97	88.8	
M1141-18	72-12018	26	1.13	9.98	9.1	90.1	
M1141-19	RBR-200055	27	1.13	9.00	8.39	92.2	
M1141-20	T872-T871	28	1.17	9.97	9.26	91.9	
M1141-22	ETGI-354	29	1.13	9.8	8.2	81.5	



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 1/15/2021

OVENTEMP IN Celsius(°C): 107
Time IN: 17:15
In Date: 01/14/2021
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:13
Out Date: 01/15/2021
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB112785

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
M1142-01	30936	9	1.00	2.00	2.00	100.0	OIL SAMPLE
M1146-01	520-Z	34	1.13	9.81	8.33	82.9	
M1146-01D	520-ZDUP	35	1.19	9.83	8.35	82.9	
M1152-01	FK-EF-01-017-ABCD	30	1.17	9.79	9.09	91.9	
M1152-02	FK-EF-01-017-ABCD	31	1.1	9.72	9.05	92.2	
M1153-01	FK-PS-01-068	32	1.13	9.71	8.27	83.2	
M1153-02	FK-PS-01-068	33	1.19	9.61	8.11	82.2	
M1156-01	VCP-EXC-SOIL	36	1.15	9.92	6.98	66.5	

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$